


Compatibilità	Codice Tecnohealth	Codice Prod.	Dimensioni	Confezione (numero pezzi)	CONNETTORE
bracciale monouso, 1 tubo, connettore maschio SLIP LUER  Compatibile: Philips HP Agilent (neonatal cuffs) ref-M18..A ; Spacelabs ref 016-....-00 , 714-....-... ; Trimline ref 28..S ; Vital Sign ref 200. , ref BP..10-20 ; Welch Allyn ref 5082-10. - . ; Zefon International ref ZDC.... ; Medline ref MDS97.. , MDS93.. ; Mabis ref MA-06-... ; Lorin Medical L00... ; Kimberly Clark ref 998.. , ref 3919. ; Invivo Research ref 9038S. ; Frontier Medical ref 15.. - . ; Ethox ref 510. ; Datascope ( Burke & Burke ) ref 0683-0. -...-... ; Critikon ( Ge Medical ) ref 21.. ; Colin ref C04..... , C040-... ; Cas ref C.. ; Baxter Allegiance Cardinal ref 30504-...S ; BMS ref 2141. ; Alaris Medical ref 503.. ;	NB-0101ML1001	BP-NN-1	neonatale, 3-6cm	10	 <p>www.tecnohealth.it                      ref NB-0101ML1XXX connector                      for all sizes</p> <p>Slip male Luer connector</p>
	NB-0101ML1002	BP-NN-2	neonatale, 4-8cm	10	
	NB-0101ML1003	BP-NN-3	neonatale, 6-11cm	10	
	NB-0101ML1004	BP-NN-4	neonatale, 7-14cm	10	
	NB-0101ML1005	BP-NN-5	neonatale, 8-15cm	10	
	bracciale monouso - 2 tubi connettore Maschio Slip Luer  Compatibile: Baxter Allegiance Cardinal ref 30504-... ; BMS ref 2142. , CAS ref 460. , C04000... , C040-0.. ; Critikon ( Ge Medical ) ref 252. , Ref 262. ; De Royal ref 7231... , ref 7231..-S ; Ethox ref	NB-0101ML2001	BP-NN2T-1	neonatale, 3-6cm	
NB-0101ML2002		BP-NN2T-2	neonatale, 4-8cm	10	

520. ; Frontier Medical ref 150. , ref 150.-S , Invivo Research ref 9038T. ; Kimberly Clark ref 3911. ; Mabis ref MA-0627. , Ge Medical Marquette ref 900005-00. ; Medline ref MDS933. , MDS975. ; Spacelabs ref 685-003.-00 , ref 714-30.-01 ; Trimline ref 280. ; Vital Sign ref 200.DT, ref BP..20-20 ; Welch Allyn ref 5082-...-2 ; Zefon International ref ZDC8.00	NB-0101ML2003	BP-NN2T-3	neonatale, 6-11cm	10
	NB-0101ML2004	BP-NN2T-4	neonatale, 7-14cm	10
	NB-0101ML2005	BP-NN2T-5	neonatale, 8-15cm	10
<b>bracciale monouso, 2 tubi / 2 tubes connettore a vite (filettato) / screw (threaded) connectors</b>  compatibile con CRITIKON	NB-0101SC2006	BP-DMDC-0	Infant, 8-14cm	5
	NB-0101SC2007	BP-DMDC-1	Pediatric, 13-20cm	5
	NB-0101SC2008	BP-DMDC-3	Small Adult, 18-26cm	5
	NB-0101SC2009	BP-DMDC-5	Adult Extra Long,26-35cm	5
	NB-0101SC2010	BP-DMDC-6	Adult, 29-38cm	5



www.tecnohealth.it  
 ref NB-0101ML2XXX  
 all sizes  
 2 tubes disposable NIBP cuffs  
 2 male slip Luer connectors



	NB-0101SC2011	BP-DMDC-7	Large Adult extra-long, 35-44cm	5
	NB-0101SC2012	BP-DMDC-8	Large Adult , 32 - 42 cm	5
	NB-0101SC2013	BP-DMDC-9	Extra Large Adult , 42-50cm	5
<b>bracciale riutilizzabile, 2 tubi                  connettore a vite maschio (male screw)</b> Compatibile: Criticare:507; Critikon:Dinamap; Invivo Research:tutti; MDE:tutti; Nihon Kohden:Lifescop L ; Spacelabs:tutti (con 2 tubi). ADC cuff ref 840-7,84, 840/1/2/3/4/5/6/6C/6I/6T/6X/; Baxter Allegiance Cardinal cuff ref 1869 , 1880 , 30432- , 30502- ,30505- ; Colin cuff ref. C04300002A- ,C043- , Critikon Ge Medical cuff ref. 27 ; De Royal cuff ref 72-10 ; Invivo Research cuff ref 9060,9070,9080,9090; Ge Medical Marquette cuff ref 9461- 301/2/3/4/5/6; Omni Cuff ref WH36 , Omron Healthcare cuff ref 21-023,21-031, Welch Allyn cuff ref 5082- -4;	NB-02020814135	BP-DMRC-0	Infant, 8-14cm	5
	NB-02021320135	BP-DMRC-1	Pediatric, 13-20cm	5
	NB-02021826135	BP-DMRC-3	Small Adult, 18-26cm	5
	NB-02022635135	BP-DMRC-5	Adult Extra Long,26-35cm	5
	NB-02022938135	BP-DMRC-6	Adult, 29-38cm	5







NIBP reusable cuffs with 2 tubes and male screw connectors

Bracciali riutilizzabili a 2 tubi connettori maschio a vite



	NB-02023242135	BP-DMRC-7	Large Adult , 32 - 42 cm	5	
	NB-02023544135	BP-DMRC-8	Large Adult extra-long, 35-44cm	5	
	NB-02024250135	BP-DMRC-9	Extra Large Adult , 42-50cm	5	
<p><b>bracciale monouso - 1 tubo,                      Connettore luer lock femmina</b></p> <p>BCI:tutti; Fukuda: tutti serie 710..                      Colin Medical monitors Pal 3100, Pal 3110P,                      Pal 3110                      Advantage 9861, Advantage 9861CO2,                      Advantage 9871, Prodigy II 2200, Prodigy II                      2210, Prodigy II 2220, Prodigy II 2240;                      Datascope:Passport XG, Passport Expert; Nihon                      Kohden: Lifescope 9, 12 Welch Allyn:Vital                      Signs/Atlas (per unità prodotte prima di Luglio                      2000) -</p>	NB-0101SDS006	BP-FLDC-0	Infant, 8-14cm	5	 <p>www.tecnohealth.it                      ref NB-0101SDSXXX                      connector                      all sizes</p> 
	NB-0101SDS007	BP-FLDC-1	Pediatric, 13-20cm	5	
	NB-0101SDS008	BP-FLDC-3	Small Adult, 18-26cm	5	
	NB-0101SDS009	BP-FLDC-5	Adult Extra Long,26-35cm	5	
	NB-0101SDS010	BP-FLDC-6	Adult, 29-38cm	5	
	NB-0101SDS011	BP-FLDC-7	Large Adult , 32 - 42 cm	5	
	NB-0101SDS012	BP-FLDC-8	Large Adult extra-long, 35-44cm	5	
	NB-0101SDS013	BP-FLDC-9	Extra Large Adult , 42-50cm	5	

<p><b>Bracciale riutilizzabile, 1 tubo connettore femmina luer lock</b></p> <p>compatibile: BCI:tutti; Fukuda serie CUF 710...                  Colin Medical monitors Pal 3100, Pal 3110P, Pal 3110                  Advantage 9861, Advantage 861CO2, Advantage 9871,                  Prodigy II 2200, Prodigy II 2210, Prodigy II 2220, Prodigy II 2240;                  Datascope: Passport XG, Passport Expert; Nihon Kohden: Lifescope 9, 12                  Welch Allyn: Vital Signs/Atlas (per unità prodotte prima di Luglio 2000) -</p> <p><b>Disponibile anche con connettore metallo sostituire 138 finale con 138M- Per prezzo con connettore metallo: richiedere ai nostri uffici</b></p>	NB-02010814138	BP-FLRC-0	Infant, 8-14cm	5	 <p>www.tecnohealth.it                  ref 138                  NIBP Reusable Cuffs 1 Tube                  Female Luer Connector                  Compatible monitors                  Datascope, Colin, BCI, Nihon                  Kohden, Welch Allyn</p> 
	NB-02011320138	BP-FLRC-1	Pediatric, 13-20cm	5	
	NB-02011826138	BP-FLRC-3	Small Adult, 18-26cm	5	
	NB-02012635138	BP-FLRC-5	Adult Extra Long, 26-35cm	5	
	NB-02012938138	BP-FLRC-6	Adult, 29-38cm	5	
	NB-02013242138	BP-FLRC-7	Large Adult , 32 - 42 cm	5	
	NB-02013544138	BP-FLRC-8	Large Adult extra-long, 35-44cm	5	
	NB-02014250138	BP-FLRC-9	Extra Large Adult , 42-50cm	5	

<p><b>Bracciale riutilizzabile 1 tubo (Per monouso ns ref NB-0101BF1XXX)</b>                      connessione baionetta                      compatibile: Criticare:508,507E Scholar II,Datascope:Accutorr Plus, Passport 2; Philips Medical:M40401 , M157.A, M455.A, tutti; Colin Medical monitors(with rectus connector) Pal 3100, Pal 3110P, Pal 3110 Advantage 9861, Advantage 9861CO2, Advantage 9871, Prodigy II 2200, Prodigy II 2210, Prodigy II 2220, Prodigy II 2240; Siemens: tutti; Spacelabs:per utilizzo con tubo 714-0018-00 e 714-0037-00 -</p> <p><b>Disponibile anche con connettore metallo sostituire 140 finale con 140M- Per prezzo con connettore metallo: richiedere ai nostri uffici</b></p>	NB-02010814140	BP-HPRC-0	Infant, 8-14cm	5	 <p>www.tecnohealth.it                      ref 140                      NIBP Reusable cuffs 1 tube                      Bayonet Connector                      Compatible monitors                      Siemens, Philips HP, Spacelabs,                      Datascope, Criticare</p> 
	NB-02011320140	BP-HPRC-1	Pediatric, 13-20cm	5	
	NB-02011826140	BP-HPRC-3	Small Adult, 18-26cm	5	
	NB-02012635140	BP-HPRC-5	Adult Extra Long, 26-35cm	5	
	NB-02012938140	BP-HPRC-6	Adult, 29-38cm	5	
	NB-02013242140	BP-HPRC-7	Large Adult , 32 - 42 cm	5	
	NB-02013544140	BP-HPRC-8	Large Adult extra-long, 35-44cm	5	

LISTINO BRACCIALI NIBP COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici  
 il presente listino può essere soggetto a variazioni senza preavviso

	NB-02014250140	BP-HPRC-9	Extra Large Adult , 42-50cm	5	
<p><b>Bracciale MONOUSO, 1 tubo                  connessione baionetta</b></p> <p>compatibile: Criticare:508,507E                  Scholar II, Datascope:Accutorr Plus, Passport2;                  Philips Medical:M40401 , M157.A, M455.A, tutti;                  Colin Medical monitors(with rectus                  connector) Pal 3100, Pal 3110P, Pal 3110                  Advantage 9861, Advantage 9861CO2,                  Advantage 9871, Prodigy II 2200, Prodigy II                  2210, Prodigy II 2220, Prodigy II 2240;                  Siemens: tutti;                  Spacelabs:per utilizzo con tubo 714-0018-00 e                  714-0037-00 -</p>	NB-0101BF1006	BP-HPDC-0	Infant, 8-14cm	5	 <p>www.tecnohealth.it                  ref. NB-0101BF1XXX                  all sizes                  1 tube NIBP disposable cuff                  bayonet connector</p> 
	NB-0101BF1007	BP-HPDC-1	Pediatric, 13-20cm	5	
	NB-0101BF1008	BP-HPDC-3	Small Adult, 18-26cm	5	
	NB-0101BF1009	BP-HPDC-5	Adult Extra Long,26-35cm	5	
	NB-0101BF1010	BP-HPDC-6	Adult, 29-38cm	5	
	NB-0101BF1011	BP-HPDC-7	Large Adult extra-long, 35-44cm	5	
	NB-0101BF1012	BP-HPDC-8	Large Adult , 32 - 42 cm	5	
	NB-0101BF1013	BP-HPDC-9	Extra Large Adult , 42-50cm	5	
		NB-02020814136	BP-MQ2TRC-0	Infant, 8-14cm	
	NB-02021320136	BP-MQ2TRC-1	Pediatric, 13-20cm	5	



<p><b>bracciali riutilizzabili a 2 tubi - utilizzabili anche con ns prolunghe ref 126 - 134/2</b></p> <p>compatibili:                      GE Medical Systems(MARQUETTE):                      compatible Ge Medical Marquette cuffs ref 00220.-E220.J , 248. - E248.J , 249. - E249.J                      compatible monitors Eagle 4000, Eagle 3000, DashTram 100, Tram 200, Tram 300. Tram 400, Tram 800,</p>	NB-02021826136	BP-MQ2TRC-3	Small Adult, 18-26cm	5	<p>ref 136                      NIBP reusable cuffs 2 tubes                      compatibles Ge Medical                      Marquette Monitors                      Eagle, Dash, Solar, Tram                      Systems</p> 
	NB-02022635136	BP-MQ2TRC-5	Adult, 26-35cm	5	
	NB-02022938136	BP-MQ2TRC-6	Adult, 29-38cm	5	
	NB-02023242136	BP-MQ2TRC-7	Large Adult, 32-42cm	5	
	NB-02023544136 -	BP-MQ2TRC-8	Large Adult extra-long, 35-44cm	5	
	NB-02024250136	BP-MQ2TRC-9	Extra Large Adult , 42-50cm	5	
<p><b>bracciali MONOUSO a 2 tubi - connessione QUICK COUPLING</b></p>	NB-0101QC2006	BP-MQ2TDC-0	Infant, 8-14cm	5	 <p>www.tecnohealth.it                      ref NB-0101QC2XXX</p>
	NB-0101QC2007	BP-MQ2TDC-1	Pediatric, 13-20cm	5	
	NB-0101QC2008	BP-MQ2TDC-3	Small Adult, 18-26cm	5	



connectors  
 all sizes

compatibili: GE Medical Systems(MARQUETTE): compatible Ge Medical Marquette cuffs ref 00220.-E220.J , 248. - E248.J , 249. - E249.J compatible monitors Eagle 4000, Eagle 3000, DashTram100, Tram 200, Tram 300. Tram 400, Tram 800,	NB-0101QC2009	BP-MQ2TDC-5	Adult Extra Long, 26-35cm	5
	NB-0101QC2010	BP-MQ2TDC-6	Adult, 29-38cm	5
	NB-0101QC2011	BP-MQ2TDC-7	Large Adult extra-long, 35-44cm	5
	NB-0101QC2012	BP-MQ2TDC-8	Large Adult , 32 - 42 cm	5
	NB-0101QC2013	BP-MQ2TDC-9	Extra Large Adult , 42-50cm	5
<b>Bracciale con sistema di gonfiaggio BULBO                      monouso, per sfigmomanometri</b>	NB-01010814BB	BP-MBDC-0	Infant, 8 – 14 cm	5
	NB-01011320BB	BP-MBDC-1	Child, 13 – 20 cm	5
	NB-01011826BB	BP-MBDC-3	Small Adult, 18 – 26 cm	5
	NB-01012635BB	BP-MBDC-5	Adult, 26 – 35 cm	5
	NB-01012938BB	BP-MBDC-6	Adult Long, 29 – 38 cm	5
	NB-010113242BB	BP-MBDC-7	Large Adult, 32-42cm	5
	NB-01013544BB	BP-MBDC-8	Large Adult Long, 35 – 44 cm	5
	NB-01014250BB	BP-MBDC-9	Extra Large Adult , 42-50cm	5
<b>Bracciale Con sistema gonfiaggio BULBO                      ritualizzabile , per sfigmomanometri</b>	NB-02020814BB	BP-MBRC-0	Infant, 8 – 14 cm	5


Senza Connettore



Senza Connettore	NB-02021320BB	BP-MBRC-1	Child, 13 – 20 cm	5
	NB-02021826BB	BP-MBRC-3	Small Adult, 18 – 26 cm	5
	NB-02022635BB	BP-MBRC-5	Adult, 26 – 35 cm	5
	NB-02023242BB	BP-MBRC-7	Large Adult, 32-42cm	5
	NB-02024250BB	BP-MBRC-9	Thigh, 42 – 50 cm	5
	NB-01022938BB	BP-MBRC-6	Adult Long, 29 – 38 cm	5
<b>Bracciale MONOUSO 1 tubo with 1 subminiature connector (quick coupler)</b>  compatibile: Welch Allyn:Protocol Propaq Ge Medical Marquette Dinamap 7000, compatibile Ge Medical Marquette Dinamap cuffs ref E9007 KF, E9007 KG, E9007KH, E9007KJ, E9007KK	NB-0101QC1006	BP-MQDC-0	Infant, 8-14cm	5
	NB-0101QC1007	BP-MQDC-1	Pediatric, 13-20cm	5
	NB-0101QC1008	BP-MQDC-3	Small Adult, 18-26cm	5
	NB-0101QC1009	BP-MQDC-5	Adult Extra Long, 26-35cm	5
	NB-0101QC1010	BP-MQDC-6	Adult, 29-38cm	5
	NB-0101QC1011	BP-MQDC-7	Large Adult extra-long, 35-44cm	5
	NB-0101QC1012	BP-MQDC-8	Large Adult , 32 - 42 cm	5



	NB-0101QC1013	BP-MQDC-9	Extra Large Adult , 42-50cm	5
<b>Bracciale riutilizzabile 1 tubo                  with 1 subminiature connector (quick coupler)</b>  compatibile: Welch Allyn:Protocol Propaq Ge Medical Marquette Dinamap 7000, compatibile Ge Medical Marquette Dinamap cuffs ref E9007 KF, E9007 KG, E9007KH, E9007KJ, E9007KK	NB-02010814142	BP-MQRC-0	Infant, 8-14cm	5
	NB-02011320142	BP-MQRC-1	Pediatric, 13-20cm	5
	NB-02011826142	BP-MQRC-3	Small Adult, 18-26cm	5
	NB-02012635142	BP-MQRC-5	Adult Extra Long,26-35cm	5
	NB-02012938142	BP-MQRC-6	Adult Long, 29 – 38 cm	5
	NB-02013242142	BP-MQRC-7	Large Adult , 32 - 42 cm	5
	NB-02013544142	BP-MQRC-8	Large Adult Long, 35 – 44 cm	5
	NB-02014250142	BP-MQRC-9	Extra Large Adult , 42- 50cm	5
	<b>bracciali riutilizzabili a 2 tubi - utilizzabili                  anche con prolunghe ns ref 128, 128/1</b> compatibili: Datex-Ohmeda:monitors tutti compatibile Ge Medical Marquette cuffs ref. 275. - E275.J and Datex Ohmeda ref cuffs 00275.. - 275..	NB-02020814137		Infant, 8-14cm
NB-02021320137			Pediatric, 13-20cm	5
NB-02021826137			Small Adult, 18-26cm	5



www.tecnohealth.it  
 ref 137  
 NIBP reusable Cuffs 2 tubes  
 Body & Insert Coupler  
 Connectors  
 Compatibles Datex Ohmeda Moni-

	NB-02022635137		Adult Extra Long, 26-35cm	5
	NB-02022938137		Adult Long, 29 – 38 cm	5
	NB-02023242137		Large Adult , 32 - 42 cm	5
	NB-02023544137		Large Adult Long, 35 – 44 cm	5
	NB-02024250137		Extra Large Adult , 42-50cm	5
bracciali riutilizzabili a 1 tubo - compatibili: Nihon Kohden :monitors tutti	NB-02010814141		Infant, 8-14cm	5
	NB-02011320141		Pediatric, 13-20cm	5
	NB-02011826141		Small Adult, 18-26cm	5
	NB-02012635141		Adult Extra Long, 26-35cm	5
	NB-02012938141		Adult Long, 29 – 38 cm	5
	NB-02013242141		Large Adult , 32 - 42 cm	5
	NB-02013544141		Large Adult Long, 35 – 44 cm	5

tors



NB-02014250141		Extra Large Adult , 42-50cm	5	
----------------	--	--------------------------------	---	--