



## Multipad

01/02

03/04

05/06

Foot

Cover



**Multipad 01**  
base area 17 x 7 cm  
variable thickness of 2 to 5 cm

**Multipad 02**  
base area 17 x 7 cm  
variable thickness of 2 to 3 cm

**Multipad Duo 01/02**  
Two connected Multipads (01 or 02)  
with one hand pump



**Multipad 03**  
base area 17 x 11 cm  
variable thickness of 2 to 5 cm

**Multipad 04**  
base area 17 x 11 cm  
variable thickness of 2 to 3 cm

**Multipad Duo 03/04**  
Two connected Multipads (03 or 04)  
with one hand pump



**Multipad 05**  
base area 33 x 28 cm  
variable thickness of 5 to 11 cm

**Multipad 06**  
base area 33 x 28 cm  
variable thickness of 5 to 7 cm

**Multipad Duo 05/06**  
Two connected Multipads (05 or 06)  
with one hand pump



**Multipad Foot 01**  
are two connected Multipads 06  
combined with a lower leg support  
which form a stable fixation for the  
foot in any position during CT or  
MR examinations.



**Cover Multipad**  
is used as a physical barrier to  
prevent a direct contact between  
the patient and the patient position-  
ing system for hygiene and comfort  
reasons. It is a single-use  
product.

## Applications



Multipad Duo 01 for CT head holder



Multipad 04 during hand assessment



Multipad Duo 06 positions the elbow



Multipad Foot 01 fixates the foot



Cover Multipad 01 applied to Multipad 01



Multipad Duo 01 for shoulder examination



Multipad Duo 03 for foot coil



Multipad Duo 06 for hand fixation



Multipad Foot 01 for foot examination



Cover applied during head MR examination



## Specific Products

### Crania

### Genus

### Carpal

### Radius

### Tibia



Crania Adult 01  
Crania Child 01

support your head examinations with perfect comfort and hygiene for the patient and lead to better image quality results, specially during time consuming assessments.



Genus Adult 01  
Genus Child 01

positioning aid is specifically designed for your most common knee scans using the standard coil. The product fits perfectly into your coil and exists in two sizes for adults and children.



Carpal Medium 01  
Carpal Large 01

are both used for forearm examinations. They comfortably position the patient's forearm and fingers and allow high-resolution image acquisition.



Radius Adult (L/R) 01

is a specific forearm positioning product designed for the XtremeCT, a high-resolution pQCT system. Its shape shows exemplarily the wide range of positioning products Pearltect is able to produce for you.



Tibia Adult 01

is a specific product designed for the XtremeCT, a high-resolution pQCT. It allows the proper fixation of the lower leg.

## Applications



Crania Adult 01 for head coil



Genus Adult 01 for knee coil



Carpal Medium 01 applied to forearm



Radius Adult 01 applied to forearm



Tibia Adult 01 applied to lower leg

