



Cavi ECG (modelli a 3 vie o 5 vie) completi di terminali pressofusi







Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	Codice Tecnohealth	STANDARD	LEADS A=US/I=EU	Defibrillation withstand	TERMINAL
AirShield	compatibile: (AIR) :HRRM-71-1/2,Porta-Fib I,II,III, System5(new),System6(Neonatal)		KV-M02010X001	AHA	3LD:A	No Resistance	SNAP (bottone)
ATL	800,800Plus,CX ATL Ultrasound CX200,CX400,Apogee,HDI 9/2000, Interspec Cardiology XL, Mark IV,Mark VII, Ultramark 4/5/7/9		KV-M02010X002	AHA	3LD:A	No Resistance	CLIP (pinzetta)
BCI	3101,3404,6100,6200, 9100,9200,Advisor, Autocorr		KV-M02010X003	IEC	3LD:I	No Resistance	SNAP (bottone)
COLIN	BP308, BP408, CBM3000		KV-M02010X004	IEC	3LD:I	No Resistance	CLIP (pinzetta)
CSI	504US,602-4USP,Poet2,Scholar		KV-M02020X001	AHA	5LD:A	No Resistance	SNAP (bottone)
GE Medical	Corometrics 556, Critikon 8100T/Dinamap 8700/9700		KV-M02020X002	AHA	5LD:A	No Resistance	CLIP (pinzetta)
IVY	101 (before 1989),402, 700		KV-M02020X003	IEC	5LD:I	No Resistance	SNAP (bottone)
LSI	LS5, 261		KV-M0202067/20	IEC	5LD:I	No Resistance	CLIP (pinzetta)
MDE	Escort 100, 101,102,200, 300, Prism, 360SL, 450SL,521A, PortaPak 80,PortaPak 90/1000		KV-M02010X005	AHA 1 KΩ	3LD:A	1KΩ	SNAP (bottone)
Medtronic / Physio Control	Lifepak 6 , VSM-1ESF (w/resistors)		KV-M02010X006	AHA 1 KΩ	3LD:A	1KΩ	CLIP (pinzetta)
Nellcor	N-200/250 C-LOCK N-1000,NPB-3900,NPB-4000		KV-M02010X007	IEC 1 KΩ	3LD:I	1KΩ	SNAP (bottone)
Nihon Kohden	B5M-8301A/8302A/8500A/8502A, OEC-6102A (s/n>20156)/7102A/8108A, OEC-6105A (s/n>20226),TEC-7100A, TEC-5200A(s/n>20421),TEC-7200A/7300A;		KV-M0201067/42-1	IEC 1 KΩ	3LD:I	1KΩ	CLIP (pinzetta)
Pace Tech	Minipack 300, Minipack 3000/3100,MiniMax4000CL, Vitalmax 800 Plus/2100/2200/4000/4100, Vitalmax/4000CL/4000 Modular/4100CL,		KV-M02020X005	AHA 1 KΩ	5LD:A	1KΩ	CLIP (pinzetta)
Philips	M3921A A1,M3923A A1, M3925AA1,M3927AA3,M3928A3, M3929AA3		KV-M02020X006	IEC 1 KΩ	5LD:I	1KΩ	CLIP (pinzetta)
Spacelabs	Burdick: Medic3/4/5/4045Defib, 500/600/900 Series (w/o respiration), 90300/90400/90600/90700 Series, TEK-208/408/511/512/514/521, PC Bedside, PCMS	KV-M02010X009	AHA	3LD:A	No Resistance	SNAP (bottone)	
TekTronix	Vitatek: 4045Defib/208/400/500Series,	KV-M02010X010	IEC	3LD:A	No Resistance	SNAP (bottone)	
Welch Allyn	Protocol: Propaq100/200Series,Atlas, Propaq Encore, Propaq CS,Ultra Smart,	KV-M02020X007	AHA	5LD:A	No Resistance	SNAP (bottone)	
Zoll	D900, PD1200, PD1400. M series	KV-M02010X011	AHA	3LD:A	No Resistance	SNAP (bottone)	
Goldway	UT4000A, UT4000A, UT4000Apro,UT6000A, UT4000B, UT4000C, UT4000F, UT6000F, UT4000Fpro, UT4800	KV-M02010X012	AHA	3LD:A	No Resistance	CLIP (pinzetta)	
Biolight Guoteng	UT4000A, UT4000A, UT4000Apro,UT6000A, UT4000B, UT4000C, UT4000F, UT6000F, UT4000Fpro, UT4800	KV-M02010X013	IEC	3LD:I	No Resistance	SNAP (bottone)	
		KV-M02010X014	IEC	3LD:I	No Resistance	CLIP (pinzetta)	

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI









 EA002S5A  EA020S5I		KV-M02020X008	AHA	5LD:A	No Resistance	SNAP (bottone)
		KV-M02020X009	AHA	5LD:A	No Resistance	CLIP (pinzetta)
		KV-M02020X010	IEC	5LD:I	No Resistance	SNAP (bottone)
		KV-M02020X011	IEC	5LD:I	No Resistance	CLIP (pinzetta)
		KV-M02010X015	AHA 1 KΩ	3LD:A	1KΩ	SNAP (bottone)
		KV-M02010X016	IEC 1 KΩ	3LD:I	1KΩ	SNAP (bottone)
		KV-M02010X017	IEC 1 KΩ	3LD:I	1KΩ	CLIP (pinzetta)
		KV-M02010X018	AHA 1 KΩ	3LD:A	1KΩ	SNAP (bottone)
		KV-M02020X012	IEC 1 KΩ	5LD:I	1KΩ	SNAP (bottone)
		KV-M0202067/43	IEC 1 KΩ	5LD:I	1KΩ	CLIP (pinzetta)
		KV-M02010X019	AHA	3LD:A	No Resistance	SNAP (bottone)
		KV-M02010X020	IEC	3LD:A	No Resistance	SNAP (bottone)
		KV-M02020X014	AHA	5LD:A	No Resistance	SNAP (bottone)
		KV-M02020X015	IEC	5LD:A	No Resistance	SNAP (bottone)

BIONET KOREA	compatibile BM3		KV-M02020X016	AHA 1KΩ	5LD:A	1KΩ	SNAP (bottone)
			KV-M02010X021	AHA 1KΩ	3LD:A	1KΩ	SNAP (bottone)
			KV-M02010X022	AHA 1KΩ	3LD:A	1KΩ	CLIP (pinzetta)
			KV-M02020X017	IEC 1KΩ	5LD:I	1KΩ	CLIP (pinzetta)
			KV-M02010X023	IEC 1KΩ	3LD:I	1KΩ	SNAP (bottone)
BIOSYS KOREA	compatibile BPM-103		KV-M02010X024	AHA 1KΩ	3LD:A	1KΩ	SNAP (bottone)
			KV-M02020X018	IEC 1KΩ	5LD:I	1KΩ	SNAP (bottone)
			KV-M02010X025	IEC 1KΩ	3LD:I	1KΩ	SNAP (bottone)
BRUKER (ODAM)	compatibile SM784, SM785 connettore con 12-pin input		KV-M02010X026	IEC	3LD:IH	No Resistance	CLIP (pinzetta)
COLIN	compatibile BP88 / BP306		KV-M02020X019	AHA	5LD:A	No Resistance	SNAP (bottone)
			KV-M02020X020	AHA	5LD:A	No Resistance	CLIP (pinzetta)
			KV-M02020X021	AHA	5LD:A	No Resistance	SNAP (bottone)
			KV-M02020X022	IEC	5LD:I	No Resistance	SNAP (bottone)
			KV-M02020X023	IEC	5LD:I	No Resistance	SNAP (bottone)
			KV-M0201067/45	IEC	3LD:I	No Resistance	CLIP (pinzetta)
			KV-M0202067/43B	IEC	5LD:I	No Resistance	CLIP (pinzetta)

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI

	compatibile BP885		KV-M02010X027	AHA	3LD:A	No Resistance	SNAP (bottone)
			KV-M02020X024	IEC	5LD:I	No Resistance	SNAP (bottone)
			KV-M02020X025	IEC	5LD:I	No Resistance	CLIP (pinzetta)
			KV-M02010X028	IEC	3LD:I	No Resistance	SNAP (bottone)
			KV-M02010X029	IEC	3LD:I	No Resistance	CLIP (pinzetta)
Creative科瑞康 Spurt世博达 Biolight宝莱特	compatibile Creative:PC9000S,PC9000,UP6000,Biolight:M12		KV-M02020X026	AHA	5LD:A	No Resistance	SNAP (bottone)
			KV-M02020X027	AHA	5LD:A	No Resistance	CLIP (pinzetta)
			KV-M02010X030	AHA	3LD:A	No Resistance	SNAP (bottone)
			KV-M02010X031	AHA	3LD:A	No Resistance	CLIP (pinzetta)
			KV-M02020X028	IEC	5LD:I	No Resistance	SNAP (bottone)
			KV-M020267/47	IEC	5LD:I	No Resistance	CLIP (pinzetta)
			KV-M020167/46	IEC	3LD:I	No Resistance	CLIP (pinzetta)
Corpuls	compatibile Corpuls 08/16, 0816M		NON DISPONIBILE	IEC	4LD:IH	No Resistance	CLIP,
DATASCOPE	2000,2001,2002,2100,2200,3000,3000A ,Ballon Pump,M/D3 w/lead select,M/D3A,M/D4, Passport,Passport IL,Passport XG,Passport XL, Passport 2,Passport 2LT,Spectrum System 90, System 90T , System 96, System 97		KV-M02020X029	AHA	5LD:A	No Resistance	SNAP (bottone)
			KV-M02020X030	AHA	5LD:A	No Resistance	CLIP (pinzetta)
			KV-M02010X032	AHA	3LD:A	No Resistance	SNAP (bottone)
			KV-M02010X033	AHA	3LD:A	No Resistance	CLIP (pinzetta)
			KV-M02020X031	IEC	5LD:I	No Resistance	SNAP (bottone)
			KV-M0202067/33	IEC	5LD:I	No Resistance	CLIP (pinzetta)
			KV-M02010X034	IEC	3LD:I	No Resistance	SNAP (bottone)
			KV-M0201067/32	IEC	3LD:I	No Resistance	CLIP (pinzetta)
DATEX	compatibile AS/3, CS/3, S/5, Cardiacap 1, Cardiacap 2 ,Cardiacap 5,Light Monitor, Satliteplus		KV-M02020X032	AHA	5LD:A	No Resistance	SNAP (bottone)
			KV-M02020X033	AHA	5LD:A	No Resistance	CLIP (pinzetta)
			KV-M02010X035	AHA	3LD:A	No Resistance	SNAP (bottone)
			KV-M02010X036	AHA	3LD:A	No Resistance	CLIP (pinzetta)
			KV-M02020X034	IEC	5LD:I	No Resistance	SNAP (bottone)
			KV-M0202067/49	IEC	5LD:I	No Resistance	CLIP (pinzetta)
			KV-M02010X037	IEC	3LD:I	No Resistance	SNAP (bottone)
			KV-M0201067/48	IEC	3LD:I	No Resistance	CLIP (pinzetta)
DRAGER SIEMEN	compatibile Drager UM1/3/3.1, RM1, PM8010, PM8014, Babyguard 8000, Cato, Babylog8000/8000 plus, PB8800 Serie 8000, Cicero, Parameterbox		KV-M0201067/22	IEC	3LD:I	No Resistance	CLIP (pinzetta)
FUKUDA DENSHI	DATASCOPE:Expert FUKUDA DENSHI:DS5100E, DS5300W, DS5400.		KV-M02010X038	AHA 1KΩ	3LD:A	1KΩ	SNAP (bottone)
			KV-M02010X039	AHA 1KΩ	3LD:A	1KΩ	CLIP (pinzetta)
			KV-M02020X035	IEC 1KΩ	5LD:I	1KΩ	SNAP (bottone)
			KV-M0202067/26	IEC 1KΩ	5LD:I	1KΩ	CLIP (pinzetta)
			KV-M02010X040	IEC 1KΩ	3LD:I	1KΩ	SNAP (bottone)
			KV-M0201067/25	IEC 1KΩ	3LD:I	1KΩ	CLIP (pinzetta)
PETAS	compatibile PETAS		KV-M02010X041	AHA	3LD:H	No Resistance	SNAP (bottone)
			KV-M02010X042	AHA	3LD:H	No Resistance	CLIP (pinzetta)
			KV-M02010X043	IEC	3LD:IH	No Resistance	SNAP (bottone)
			KV-M02010X044	IEC	3LD:IH	No Resistance	CLIP (pinzetta)

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso










LISTINO PREZZI RICAMBI COMPATIBILI

GE-MEDICAL MARQUETTE	Compatibile GEEMS:Dash PRO 4000,Dash PRO 3000, Dash PRO 2000,Dash 1000, Eagle Systems, MAC-Lab, Solar 9500, Solar 8000. Solar 3000, Tram 100, Tram 200,Tram 300, Tram 400, Tram 500,Tram 600, Tram 800, Unity Systems Transport Pro, Ecografi VIVID e LOGIQ series		KV-M02020X036	AHA	5LD:A	No Resistance	SNAP (bottone)
			KV-M02020X037	AHA	5LD:A	No Resistance	CLIP (pinzetta)
			KV-M02010X045	AHA	3LD:A	No Resistance	SNAP (bottone)
			KV-M02010X046	AHA	3LD:A	No Resistance	CLIP (pinzetta)
			KV-M02020X038	IEC	5LD:I	No Resistance	SNAP (bottone)
			KV-M0202067/28	IEC	5LD:I	No Resistance	CLIP (pinzetta)
			KV-M02010X047	IEC	3LD:I	No Resistance	SNAP (bottone)
			KV-M0201067/27	IEC	3LD:I	No Resistance	CLIP (pinzetta)
			KV-M02010X048	AHA 1KΩ	3LD:A	1KΩ	SNAP (bottone)
			KV-M02010X049	IEC 1KΩ	3LD:I	1KΩ	SNAP (bottone)
	compatibile PRO1000		KV-M02010X050	IEC 1KΩ	3LD:I	1KΩ	CLIP (pinzetta)
GE-MEDICAL	compatibile GE MEDICAL SYSTEMS:2000, 2020, 2110, 2110N, 2120, 2210 2230, 2230N, 2310, 2330, 2330N, 7000 Series, Tram AR		KV-M02020X039	AHA	5LD:A	No Resistance	SNAP (bottone)
			KV-M02020X040	AHA	5LD:A	No Resistance	CLIP (pinzetta)
			KV-M02010X051	AHA	3LD:A	No Resistance	SNAP (bottone)
			KV-M02010X052	AHA	3LD:A	No Resistance	CLIP (pinzetta)
			KV-M02020X041	IEC	5LD:I	No Resistance	SNAP (bottone)
			KV-M02020X042	IEC	5LD:I	No Resistance	CLIP (pinzetta)
			KV-M02010X053	IEC	3LD:I	No Resistance	SNAP (bottone)
			KV-M02010X054	IEC	3LD:I	No Resistance	CLIP (pinzetta)
			KV-M02020X043	AHA	5LD:A	No Resistance	SNAP (bottone)
			KV-M02020X044	AHA	5LD:A	No Resistance	CLIP (pinzetta)
HELLIGE	compatibile Hellige,Servomed 4 , Generation, SCP 851, SCP 852 , SCP 844,SCS 804, CardioServ , VICOM-sm		KV-M02010X055	AHA	3LD:A	No Resistance	SNAP (bottone)
			KV-M0202067/22d	IEC	5LD:I	No Resistance	CLIP (pinzetta)
			KV-M0201067/22c	IEC	3LD:I	No Resistance	CLIP (pinzetta)
			KV-M02020X045	AHA	5LD:A	No Resistance	SNAP (bottone)
KONTRON	compatibile 7000, 7135, 7135B, 7141, 7210 7250, 7251, K2000, KAAAT, Kolormon, Micromon, Minimon, Supermon , Vitalmon		KV-M02020X046	AHA	5LD:A	No Resistance	CLIP (pinzetta)
			KV-M02020X047	IEC	5LD:I	No Resistance	SNAP (bottone)
			KV-M0202067/44a	IEC	5LD:I	No Resistance	CLIP (pinzetta)
			KV-02010332	IEC	3LD:I	No Resistance	SNAP (bottone)
			KV-M0201067/44	IEC	3LD:I	No Resistance	CLIP (pinzetta)
			KV-M02010X056	AHA	3LD:A	No Resistance	CLIP (pinzetta)
			KV-M02010X057	IEC	3LD:I	No Resistance	SNAP (bottone)
			KV-M02010X058	IEC	3LD:I	No Resistance	CLIP (pinzetta)
			KV-M02010X059	AHA 1KΩ	3LD:H	1KΩ	SNAP (bottone)
			KV-M02010X060	AHA 1KΩ	3LD:H	1KΩ	CLIP (pinzetta)
L&T (INDIA)	compatibilità Dlanet50		KV-M02010X061	IEC 1KΩ	3LD:IH	1KΩ	SNAP (bottone)
			KV-M02010X062	IEC 1KΩ	3LD:IH	1KΩ	CLIP (pinzetta)
			KV-M02010X063	AHA	3LD:A	No Resistance	SNAP (bottone)
			KV-M02010X064	AHA	3LD:A	No Resistance	CLIP (pinzetta)
	compatibilità MP1000,MP600,MP500		KV-M02020X049	IEC	5LD:I	No Resistance	CLIP (pinzetta)
			KV-M02010X065	IEC	3LD:I	No Resistance	SNAP (bottone)
			KV-M02010X066	IEC	3LD:I	No Resistance	CLIP (pinzetta)
			KV-M02010X067	AHA	3LD:A	No Resistance	SNAP (bottone)
compatibilità M160246 M360697		KV-M02010X068	AHA	3LD:A	No Resistance	CLIP (pinzetta)	
		KV-M02020X050	IEC	5LD:I	No Resistance	CLIP (pinzetta)	
		KV-M02010X069	IEC	3LD:I	No Resistance	CLIP (pinzetta)	

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI

M&B			KV-M02010X070	AHA	3LD:A	No Resistance	SNAP (bottone)
			KV-M02020X051	AHA	5LD:A	No Resistance	SNAP (bottone)
MENNEN	compatibilità Horizon 250,3000,450, Microsolo,936 Series,425/101 Preamp,Cardio-Pak 2000,740,741,742		KV-M02010X071	AHA 1KΩ	3LD:A	1KΩ	SNAP (bottone)
			KV-M02010X072	AHA 1KΩ	3LD:A	1KΩ	CLIP (pinzetta)
			KV-M02010X073	IEC 1KΩ	3LD:I	1KΩ	SNAP (bottone)
			KV-M02010X074	IEC 1KΩ	3LD:I	1KΩ	CLIP (pinzetta)
	compatibilità (MENNEN):900,1100,1110, Horizon 2000,2100,2200,5200 Horizon 9000XL,MR1200,1270 Horizon 1280,1300,Horizon 1380		KV-M02010X075	AHA	3LD:A	No Resistance	SNAP (bottone)
			KV-M02010X076	AHA	3LD:A	No Resistance	CLIP (pinzetta)
			KV-M02020X052	IEC	5LD:I	No Resistance	SNAP (bottone)
			KV-M0202067/53	IEC	5LD:I	No Resistance	CLIP (pinzetta)
			KV-M0201067/52	IEC	3LD:I	No Resistance	CLIP (pinzetta)
			KV-M02010X077	AHA	3LD:A	No Resistance	SNAP (bottone)
MINDRAY	compatibile PM9000,PM8000,PM7000 MEC1000, MEC1200 connettore 6-pin		KV-M02010X078	AHA	3LD:A	No Resistance	CLIP (pinzetta)
			KV-M02010X079	IEC	3LD:I	No Resistance	SNAP (bottone)
			KV-M0201067/20/A	IEC	3LD:I	No Resistance	CLIP (pinzetta)
			KV-M02020X053	AHA	5LD:A	No Resistance	SNAP (bottone)
			KV-M02020X054	AHA	5LD:A	No Resistance	CLIP (pinzetta)
			KV-M02020X055	IEC	5LD:I	No Resistance	SNAP (bottone)
			KV-M0202067/20	IEC	5LD:I	No Resistance	CLIP (pinzetta)
			compatibile BeneView (T5, T6, T8) IPM series connettore (diritto) 12-pin , T8		KV-M02020X056	AHA	5LD:A
	KV-M02010X080	AHA			3LD:A	No Resistance	SNAP (bottone)
	KV-M02010X081	AHA			3LD:A	No Resistance	CLIP (pinzetta)
KV-M02020X057	IEC	5LD:I			No Resistance	SNAP (bottone)	
KV-M0202067/20B	IEC	5LD:I			No Resistance	CLIP (pinzetta)	
KV-M02010X082	IEC	3LD:I			No Resistance	SNAP (bottone)	
KV-M02020X058	AHA	5LD:A			No Resistance	SNAP (bottone)	
KV-M02010X083	AHA	3LD:A			No Resistance	SNAP (bottone)	
compatibile PM5000/6000 connettore (angolato) 12-pin , T8		KV-M02020X059	IEC	5LD:I	No Resistance	SNAP (bottone)	
		KV-M02010X084	IEC	3LD:I	No Resistance	SNAP (bottone)	
		KV-M0201067/20C	IEC	3LD:I	No Resistance	CLIP (pinzetta)	
		KV-M02020X060	AHA	5LD:A	No Resistance	SNAP (bottone)	
NEC	compatibile 3M12		KV-M02020X061	AHA	5LD:A	No Resistance	CLIP (pinzetta)
			KV-M02010X085	AHA	3LD:A	No Resistance	SNAP (bottone)
			KV-M02010X086	IEC	3LD:I	No Resistance	SNAP (bottone)
			KV-M02010X087	AHA	3LD:A	No Resistance	CLIP (pinzetta)
	compatibile PB1916E		KV-M02010X088	IEC	3LD:I	No Resistance	SNAP (bottone)
			KV-M0201067/29	IEC 1KΩ	3LD:I	1KΩ	SNAP (bottone)
NIHON KOHDEN	Compatibile Nihon Kodhen TEK con connettore circolare a 8 PIN compatibile OEC-6102A (s/n<20156), OEC-6105A (s/n<20226), OEC-6120A/J/K BSM-7105K		KV-M02010X089	IEC 1KΩ	3LD:I	1KΩ	CLIP (pinzetta)

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI

	compatibile BSM-7102A,OMP-7201C,OMP-7203C, TEC-5200A (s/n<20421) Connettore lato monitor/Defibrillatore 11 Pin diam. esterno connettore 20mm – NO resistenza -		KV-M0201067/31-1	IEC	3LD:I	No Resistance	SNAP (bottone)			
			KV-M0201067/31	IEC	3LD:I	No Resistance	CLIP (pinzetta)			
	compatibile TEC-7621C, TEC-7621C Connettore lato monitor/Defibrillatore 11 Pin diam. Connettore diam . 20mm		KV-M0201067/31B	IEC	3LD:I	No Resistance	CLIP (pinzetta)			
			KV-M02010X090	AHA 200KΩ	3LD:A	200KΩ	CLIP (pinzetta)			
	compatibile TEC-7621K, TEC-7631K Connettore diam . 21,5mm		KV-M02010X091	IEC 100KΩ	3LD:I	100KΩ	SNAP (bottone)			
			KV-M0201067/31A	IEC 100KΩ	3LD:I	100KΩ	CLIP (pinzetta)			
PETAS	compatibile KMA160R,165R,260R,K-250,300...		KV-M02010X092	AHA	3LD:H	No Resistance	SNAP (bottone)			
			KV-M02010X093	AHA	3LD:H	No Resistance	CLIP (pinzetta)			
			KV-M02010X094	IEC	3LD:IH	No Resistance	SNAP (bottone)			
			KV-M02010X095	IEC	3LD:IH	No Resistance	CLIP (pinzetta)			
PHILIPS	compatibile 43100A,43110A,43120A,78352C 78354C,78834C,862474 C3, M1001A,M1002A,M1165A,M1166A,M1167A,M1175A, M1177A,M1205A,M1264A,M1275A Transport M1280A,M1722A/B,M1723A/B M2475B,M2601A,M3000A,M3001A M3046A M2/M3/M4,M4735A connettore circolare 12 pin		KV-M02020X062	AHA 1KΩ	5LD:A	1KΩ	SNAP (bottone)			
			KV-M02020X063	AHA 1KΩ	5LD:A	1KΩ	CLIP (pinzetta)			
			KV-M02010X096	AHA 1KΩ	3LD:A	1KΩ	SNAP (bottone)			
			KV-M02010X097	AHA 1KΩ	3LD:A	1KΩ	CLIP (pinzetta)			
			KV-M02020X064	IEC 1KΩ	5LD:I	1KΩ	SNAP (bottone)			
			KV-M0202067/33B1	IEC 1KΩ	3LD:I	1KΩ	CLIP (pinzetta)			
			KV-M02010X098	IEC 1KΩ	3LD:I	1KΩ	SNAP (bottone)			
			KV-M0201067/33B	IEC 1KΩ	3LD:I	1KΩ	CLIP (pinzetta)			
			KV-M02020X065	AHA 1KΩ	5LD:A	1KΩ	SNAP (bottone)			
			KV-M02020X066	AHA 1KΩ	5LD:A	1KΩ	CLIP (pinzetta)			
			KV-M02010X099	AHA 1KΩ	3LD:A	1KΩ	SNAP (bottone)			
			KV-M02010X100	AHA 1KΩ	3LD:A	1KΩ	CLIP (pinzetta)			
			KV-M02020X067	IEC 1KΩ	5LD:I	1KΩ	SNAP (bottone)			
			KV-M0202067/33A1	IEC 1KΩ	3LD:I	1KΩ	CLIP (pinzetta)			
Primedic			KV-M02010X102	AHA 10KΩ	3LD:H	10KΩ	SNAP (bottone)			
			KV-M02010X103	AHA 10KΩ	3LD:H	10KΩ	CLIP (pinzetta)			
			KV-M02010X104	IEC 10KΩ	3LD:IH	10KΩ	SNAP (bottone)			
			KV-M02010X105	IEC 10KΩ	3LD:IH	10KΩ	CLIP (pinzetta)			
			MEDTRONIC PHYSIOCONTROL	compatibile Lifepak 12 , 120, Lifepak 20		KV-M02020X068	AHA 10KΩ	5LD:A	10KΩ	SNAP (bottone)
						KV-M02020X069	AHA 10KΩ	5LD:A	10KΩ	CLIP (pinzetta)
KV-M0202067/24						No Resistance				
KV-M02010X106	AHA 10KΩ	3LD:A				10KΩ	SNAP (bottone)			
KV-M02010X107	AHA 10KΩ	3LD:A				10KΩ	CLIP (pinzetta)			
KV-M02020X070	IEC 10KΩ	5LD:I				10KΩ	SNAP (bottone)			
MEDTRONIC PHYSIOCONTROL	compatibile PHYSIO-CONTROL:Lifepak 2,3,33. Scope Monitor II ; SIEMENS:Porta Pak 500PB, Porta Pak 500PM ; SIMONSEN & WEELS(S&W);CB-600 Defib, DMS-600 ; SURVIVAL INSTRUMENTS:Safeguard I, II UNIRAD:Series C Echo,VIGGOSPECTRAMED (Gould):Ekoline 20A, 21, Smith-Kline 20A		KV-M02010X109	AHA	3LD:A	No Resistance	SNAP (bottone)			
			KV-M02010X110	AHA	3LD:A	No Resistance	CLIP (pinzetta)			
			KV-M02010X111	IEC	3LD:I	No Resistance	SNAP (bottone)			
			KV-M02010X112	IEC	3LD:I	No Resistance	CLIP (pinzetta)			

LISTINO CAVI ECG COMPATIBILI


PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI

SAADAT			KV-M02010X113	IEC 1KΩ	3LD:IH	1KΩ	SNAP (bottone)
			KV-M02010X114	IEC 1KΩ	3LD:IH	1KΩ	SNAP (bottone)
S&W (Artema) Diascope	compatibile con: Artema S&W mod. Neoscope, Diascope mod 1 e 2, Diascope NT, Defibrillator, Athena System, DSNT, Traveller, 8000 System, DS1, DS2, DMS 700;		KV-M02010X115	AHA	3LD:A	No Resistance	SNAP (bottone)
			KV-M02020X071	IEC	5LD:I	No Resistance	CLIP (pinzetta)
			KV-M02010X116	IEC	3LD:I	No Resistance	SNAP (bottone)
			KV-M020167/21	IEC	3LD:I	No Resistance	CLIP (pinzetta)
SIEMENS	compatibile SIEMENS:400 Series, 500 Series, 600 series, 700 Series, 900 Series, 1200 Series Sirecust mod. 341, 400, 404, 500 series, 610, 620, 630, 720, 722, 730, 732, 900, 960, 961, 1260, 1281, 1801, Polymed Module, Modul Box Sirem.		KV-M0202067/41	IEC	5LD:I	No Resistance	CLIP (pinzetta)
			KV-M0201067/40	IEC	3LD:I	No Resistance	CLIP (pinzetta)
SCHILLER			KV-M02020X072	IEC	5LD:I	No Resistance	CLIP (pinzetta)
SPACELABS	compatibile 90496 Ultraview, 90369 Ultraview, 90367 Ultraview		KV-M02020X073	AHA	5LD:A	No Resistance	SNAP (bottone)
			KV-M02020X074	AHA	5LD:A	No Resistance	CLIP (pinzetta)
			KV-M02010X117	AHA	3LD:A	No Resistance	CLIP (pinzetta)
			KV-M02020X075	IEC	5LD:I	No Resistance	SNAP (bottone)
			KV-M0202067/51	IEC	5LD:I	No Resistance	CLIP (pinzetta)
			KV-M02010X118	IEC	3LD:I	No Resistance	SNAP (bottone)
			KV-M0201067/50	IEC	3LD:I	No Resistance	CLIP (pinzetta)
SPACELABS BURDICK	GE MEDICAL SYSTEMS:Corometrics 503, Corometrics 504,Corometrics 508; SIEMENS ELECTROMEDICAL SYSTEMS:806180A Nuclear Camera		KV-M02010X119	AHA	3LD:H	No Resistance	SNAP (bottone)
			KV-M02010X120	AHA	3LD:H	No Resistance	CLIP (pinzetta)
			KV-M02010X121	IEC	3LD:IH	No Resistance	SNAP (bottone)
			KV-M02010X122	IEC	3LD:IH	No Resistance	CLIP (pinzetta)
WELCH ALLYN	compatibile Propaq 100 Series, Propaq CS , Ultra Smart , Atlas,Propaq 200 Series, Propaq Encore,		KV-M02020X076	AHA	5LD:A	No Resistance	SNAP (bottone)
			KV-M02010X123	AHA	3LD:A	No Resistance	SNAP (bottone)
			KV-M02020X077	IEC	5LD:I	No Resistance	SNAP (bottone)
			KV-M0202067/43	IEC	5LD:I	No Resistance	CLIP (pinzetta)
			KV-M02010X124	IEC	3LD:I	No Resistance	CLIP (pinzetta)
			KV-M0201067/42	IEC	3LD:I	No Resistance	CLIP (pinzetta)
















LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso



Yoke type	Productor guide	TH code	Description	Compatible Brand	Need Adapter
DIN type, Code: C		KV-S02010X001	3 vie, 2 metri, AAMI Leads	Datascope Welch Allyn	EC002-3AI
		KV-S02010X002	3 vie, 2 metri, IEC Leads	Spacelabs	EC020-3AI
				Biosys	EC040-3AI
			KV-S02020X001	5 vie, 2,5 metri, AAMI Leads	Creative
	MEK				
	KV-S02020X002		5 vie, 2,5 metri, IEC Leads	Lohmeier	EC078-3AI
				Datascope	EC002-5AI
	AAMI type, Code: D		KV-S02010X003	3 vie, 2 metri, AAMI Leads	PHILIPS
PHILIPS					ED035-3AI
KV-S02010X004			3 vie, 2 metri, IEC Leads	PHILIPS	ED036-3AI
		KV-S02020X003	5 vie, 2,5 metri, AAMI Leads	PHILIPS	ED023-5AI
				PHILIPS	ED035-5AI
		KV-S02020X004	5 vie, 2,5 metri, IEC Leads	PHILIPS	ED036-5AI
VS type, Code: E		KV-S02010X005	3 vie, 2 metri, AAMI Leads	Marquette	EE007-5AI
		KV-S02010X006	3 vie, 2 metri, IEC Leads		


LISTINO PREZZI RICAMBI COMPATIBILI

		KV-S02020X005	5 vie, 2,5 metri, AAMI Leads		Marquette	EE007-5AI		
					Hellige	EE011-5AI		
		KV-S02020X006	5 vie, 2,5 metri, IEC Leads					
Euro type, Code: J		KV-S02010X007	3 vie, 2 metri, AAMI Leads		Siemens	EJ052-3AI		
					Mindray	EJ021-3AI		
		KV-S02010X008	3 vie, 2 metri, IEC Leads		Mindray	EJ050-3AI		
					Mindray	EJ051-3AI		
		KV-S02020X007	5 vie, 2,5 metri, AAMI Leads		Siemens	EJ052-5AI		
					Mindray	EJ021-5AI		
		KV-S02020X008	5 vie, 2,5 metri, IEC Leads		Mindray	EJ050-5AI		
					Drager	EJM079-5AI		
TRU-LINK type, Code: L		KV-S02010X009	3 VIE, 2 metri, AAMI Leads		Fukuda	VJ032-EAI		
					Lohmeter	EJM079-5AI		
		KV-S02010X010	3 VIE, 2 metri, IEC Leads					
		KV-S02020X009	5 VIE, 2,5 metri, AAMI Leads		Spacelabs	EL009-5AI		
				Spacelabs	EL009-5AI			
		KV-S02020X010	5 VIE, 2,5 metri, IEC Leads					

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso

ECG Trunk Cable and Leads for Patient monitor



TRUNK CABLE - CAVO BASE

Compatible Brand	Compatible Models No.	OEM P/N	TH CODE	DESCRIPTION (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	YOKE (LEADS SIDE)
Din Style Normal	Compatibile AIR SHIELDS (NARCO,Vickers): HRRM-71-1/2 ,Porta-Fib I, II, III, System 5(new),System 6 (Neonatal); ATL ULTRASOUND:800 ,800 Plus,Apogee , CX ATL ULTRASOUND: CX200,CX400,HDI 9/2000,Interspec	D-1340	KV-0201067/12	Din 3Ld Patient Cable,3 metri AHA/IEC, straight 6pin>Din 3J, no resistance fit for Din style LeadWires	004	101
		D-1540	KV-02020X001	Din 5Ld Patient Cable,3 metri AHA/IEC, straight 6pin>Din 5J, no resistance fit for Din style LeadWires	004	102
			KV-02010X001	Din 3Ld Patient Cable,3 metri AHA/IEC, straight 6pin>Din 3J, 1KΩ resistance fit for Din style LeadWires	003	101
			KV-02020X002	Din 5Ld Patient Cable,3 metri AHA/IEC, straight 6pin>Din 5J, 1KΩ resistance fit for Din style LeadWires	003	102
Lifelinc (LL) Style Normal	compatibile LIFE SCIENCE INSTRUMENTS (LSI):LS5 , 261 MEDICAL DATA electronics (MDE):Escort 100, Escort 101, Escort 102 , Escort 200,Escort 300, Prism , 360SL, 450SL, 521A PortaPak 80 , PortaPak 90 , PortaPak 1000;		KV-02010038/1	LL 3Ld Patient Cable,3 metri AHA/IEC, straight 6pin>LL 3LD, no resistance fit for LL style LeadWires	004	107
			KV-02020X003	LL 5Ld Patient Cable,3 metri AHA/IEC, straight 6pin>LL 5LD, no resistance fit for LL style LeadWires	004	105
			KV-02010X003	LL 3Ld Patient Cable,3 metri AHA/IEC, straight 6pin>LL 3LD, 1KΩ resistance fit for LL style LeadWires	003	107
			KV-02020X004	LL 5Ld Patient Cable,3 metri AHA/IEC, straight 6pin>LL 5LD, 1KΩ resistance , fit for LL style LeadWires	003	105
Conmed Style Normal	PHILIPS MEDICAL SYSTEMS: M3921A A1 Compact Monitor , M3923A A1 Compact Monitor, M3925A A1 Compact Monitor ,		KV-02010X004	Conmed 3Ld Patient Cable, Length: 3 metri ; AHA/IEC, straight 6pin>Conmed 3LD, no resistance, fit for Conmed style LeadWires	004	671

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso














LISTINO PREZZI RICAMBI COMPATIBILI

	M3927A A3 Compact Monitor , M3929A A3 Compact Monitor PIONEER (American Optical): 26125 , 33167, 5260 , 5330 ,Pulsar 4;		KV-02020X005	Conmed 5Ld Patient Cable, Length: 3 metri ; AHA/IEC, straight 6pin>Conmed 5LD, no resistance, fit for Conmed style LeadWires		
			KV-02010X005	Conmed 3Ld Patient Cable, Length: 3 metri ; AHA/IEC, straight 6pin>Conmed 3LD, 1KQ resistance, fit for Conmed style LeadWires		
			KV-02020X006	Conmed 5Ld Patient Cable, Length: 3 metri ; AHA/IEC, straight 6pin>Conmed 5LD, 1KQ resistance, fit for Conmed style LeadWires		
Bionet Korea	compatibile BM3		KV-02010X006	3Ld Patient Cable,3 metri AHA/IEC, Round 8pin>LL-3LD, 1KQ resistance fit for LL style LeadWires		
Biosys Korea	compatibile BPM-103		KV-02010X007	3Ld Patient Cable,3 metri AHA/IEC, Readel 6pin>Din 3J, 1KQ resistance fit for Din style LeadWires		
			KV-02020X007	5Ld Patient Cable,3 metri AHA/IEC, Readel 6pin>Din 5J, 1KQ resistance fit for Din style Lead		
Colin	compatibile BP88 / BP306		KV-02020X008	5Ld Patient Cable,3 metri IEC, L type round 6pin>6J, no resistance fit for Colin LeadWires		
	compatibile BP88S		KV-02020X008/1	5Ld Patient Cable,3 metri IEC, straight 6P>6J, no resistance fit for Colin LeadWires		
Creative Spurt Biolight	compatibile Creative: PC9000S,PC9000,UP6000, Biolight: M12		KV-02010X008	3Ld Patient Cable,3 metri AHA/IEC, Readel 6pin>Din 3J, 1KQ resistance fit for Din style LeadWires		
			KV-02010XX18	5Ld Patient Cable,3 metri AHA/IEC, Readel 6pin>Din 5J, 1KQ resistance fit for Din style LeadWires		
Datascope	compatibile DATASCOPE 2000,2001,2002,2100,2200,3000, 3000A ,Ballon Pump,M/D3 w/lead select,M/D3A,M/D4,Passport, Passport IL,Passport XG,Passport XL,Passport 2,Passport 2LT,Spectrum System 90, System 90T , System 96, System 97,System 98, System 98XT, Trio		KV-02010067/19	Din 3Ld Patient Cable,3 metri AHA/IEC, straight 6pin>Din 3J, no resistance fit for Din style LeadWires		
			KV-02010XX19	Din 5Ld Patient Cable,3 metri AHA/IEC, straight 6pin>Din 5J, no resistance fit for Din style LeadWires		
Datex	compatibile AS/3, CS/3, S/5, Cardiocap 1, Cardiocap 2 , Cardiocap 5 , Light Monitor , Satliteplus		KV-02010061	3Ld Patient Cable,3 metri AHA/IEC, Round 10I>LL 3LD, no resistance fit for LL style LeadWires		
			KV-02020063	5Ld Patient Cable,3 metri AHA/IEC, Round 10I>LL 5LD, no resistance fit for LL style LeadWires		

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI

Datex-Ohmeda	compatibile Datex-Ohmeda Patient monitor		KV-0201015	3Ld Patient Cable,3 metri AHA/IEC, Round 10I>ASP 3LD, 4.7KΩ resistance fit for ASP style LeadWires			183
			KV-0201017	5Ld Patient Cable,3 metri AHA/IEC, Round 10I>ASP 5LD, 4.7KΩ resistance fit for ASP style LeadWires			221
Drager Siemens	compatibile Drager UM1/3/3.1, RM1, PM8010, PM8014, Babyguard 8000, Cato, Babylog8000/8000 plus, PB8800,Serie 8000, Cicero, Parameterbox,		KV-0201067	5G 3Ld Patient Cable,3 metri AHA/IEC, Round 10Pin>Din 3I, no resistance fit for Din style LeadWire			101
	compatibile Multimed-Pot system fur serie SC 6000, SC 6002XL, SC 7000, SC 9000, Art. Nr.3368391 (8.2ft) Art. Nr. 5950196 (4.3 metri)	3368391	KV-030171/A	5C9000XL Multi-link Cable, 8.2ft AHA/IEC, 1KΩ resistance fit for Euro style LeadWires			179
Siemens	compatibile Sirecust Series 400, 500, 600, Serie 700, 900, 1200, Sirecust 341/610/620/630/720/722, Sirecust 730/732/900/960/961/1260, Sirecust 1281/1801/400/-500, ECG Monitor, Polymed Modul, Modulbox Sirem		KV-02010030/2	3Ld Patient Cable,3 metri AHA/IEC, Round 10I>Euro 3LD, no resistance fit for Euro style LeadWires			242
			KV-02010XX21	5Ld Patient Cable,3 metri AHA/IEC, Round 10I>Euro 5LD, no resistance fit for Euro style LeadWires			243
GE-Medical Marquette	compatibile GEMS:Dash PRO 4000, Dash PRO 3000, Dash PRO 2000, Dash 1000, Eagle Systems, MAC-Lab, Solar 9500, Solar 8000, Solar 3000, Tram 100, Tram 200, Tram 300, Tram 400, Tram 500, Tram 600, Tram 800, Unity Systems		KV-02020001	5Ld Patient Cable,3 metri AHA, (IEC) Round 10I>VS 5LD, no resistance fit for VS style LeadWires			089
GE-Hellige	Cardioserv Eagle 1000/3000 Hellige,Servomed 4 , Generation, SCP 851, SCP 852 , SCP 844, SCS 804,VICOM-sm		KV-02020X010	5Ld Patient Cable,3 metri AHA, (IEC) Round 10I>VS 5LD, 10KΩ resistance fit for VS style LeadWires			089
			KV-02020X011	5Ld Patient Cable,3 metri AHA/IEC, Round 10P>L 5LD, no resistance fit for LL style LeadWires			105
			KV-02010X009	3Ld Patient Cable,3 metri AHA/IEC, Round 10I>ASP 3LD, 1KΩ resistance fit for ASP style LeadWires			183I
			KV-0202067/B	5Ld Patient Cable,3 metri AHA/IEC, Round 10I>ASP 5LD, 1KΩ resistance fit for ASP style LeadWires			221
GE-Medical Marquette	compatibile Pro1000		KV-02010X010	3Ld Patient Cable,3 metri AHA/IEC, Round 6pin>ASP 3LD, 1KΩ resistance fit for ASP style LeadWires			183I
			KV-02020X012	5Ld Patient Cable,3 metri AHA/IEC, Round 6pin>ASP 5LD, 1KΩ resistance fit for ASP style LeadWires			221
Goldway	compatibile UT4000A,UT4000Apro,UT6000A, UT4000B, UT4000C, UT4000E		KV-02010X011	3Ld Patient Cable,3 metri AHA/IEC, Round 6p>LL 3LD, no resistance fit for LL style LeadWires			10r

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso





























LISTINO PREZZI RICAMBI COMPATIBILI

	UT6000F, UT4000Fpro, UT4800		KV-02020X013	5Ld Patient Cable,3 metri AHA/IEC, Round 6p>LL 5LD, no resistance fit for LL style LeadWires		
Huntleigh Healthcare	compatibile UK MONITOR		KV-02010X012	Din 3Ld Patient Cable,3 metri AHA/IEC, DB9M>Din 3I, no resistance fit for Din style LeadWires		
Kontron	compatibile KONTRON 7000, 7135, 7135B, 7141, 7210,7250, 7251, K2000, KAAT, Kolormon, Micromon, Minimon,Supermon , Vitalmon		KV-02020X014	5Ld Patient Cable,3 metri AHA/IEC, Round 12P>LL 5LD, no resistance fit for LL style LeadWires		
MEK Korea	compatibile MP1000,MP600,MP500		KV-02020X015	Din 5Ld Patient Cable,3 metri AHA/IEC, DB9M>Din 5I, no resistance fit for Din style LeadWires		
MINDRAY	compatibile PM-7000, PM-8000series, PM-9000series, new PM-6000	0010-30-12378	KV-02010X013	3Ld Patient Cable,3 metri AHA/IEC, Round 6pin>Euro 3LD, 1KΩ resistance fit for Euro style LeadWires		
		0010-30-12257	KV-02020X016	5Ld Patient Cable,3 metri AHA/IEC, Round 6pin>Euro 5LD, 1KΩ resistance fit for Euro style LeadWires		
	compatibile PM-6800,2004年以后的PM-6000, connettore (diritto) 12-pin verde, BeneView (T5, 56, T8), iPM Series	0010-30-42724	KV-02010X014	3Ld Patient Cable,3 metri AHA/IEC, straight Round 12P>Euro 3LD, 1KΩ resistance /inductance fit for Euro style LeadWires		
		0010-30-42723	KV-02020X017	5Ld Patient Cable,3 metri AHA/IEC, straight Round 12P>Euro 5LD, 1KΩ resistance /inductance fit for Euro style LeadWires		
	PM6000, connettore (diritto) 12-pin verde,T8	0010-30-12379	KV-02010X015	3Ld Patient Cable,3 metri AHA/IEC, Type L Round 12P>Euro 3LD, 1KΩ resistance fit for Euro style LeadWires		
		0010-30-12260	KV-02020X018	5Ld Patient Cable,3 metri AHA/IEC, Type L Round 12P>Euro 5LD, 1KΩ resistance fit for Euro style LeadWires		
	PM6000, connettore (angolato) 12-pin verde		KV-02010X016	3Ld Patient Cable,3 metri AHA/IEC, Type L Round 12P>Euro 3LD, 1KΩ resistance fit for Euro style LeadWires		
			KV-02020X019	5Ld Patient Cable,3 metri AHA/IEC, Type L Round 12P>Euro 5LD, 1KΩ resistance fit for Euro style LeadWires		
NEC	compatibile 3M12	47524	KV-02010X017	3Ld Patient Cable,3 metri AHA/IEC, Round 16p>NEC 3LD, no resistance fit for 47504 LeadWires		
		47525	KV-02020X020	5Ld Patient Cable,3 metri AHA/IEC, Round 16p>NEC 5LD, no resistance fit for 47505 LeadWires		

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso























LISTINO PREZZI RICAMBI COMPATIBILI

Nihon Kohden	compatibile OEC-6102A (s/n<20156), OEC-6105A (s/n<20226), OEC-6120A/J/K BSM-7105K		KV-0201067/12	3Ld Patient Cable,3 metri AHA/IEC, Round 8p>Din 3LD, no resistance fit for Din style LeadWires		237		101
		JC-005P	KV-02010X018	3Ld Patient Cable,3 metri IEC, Round 8p>NK 6J, no resistance fit for BR-019 LeadWires		237		155
Nihon Kohden	compatibile BSM-7102A,OMP-7201C, OMP-7203C, TEC-5200A (s/n<20421),		KV-02010X019	3Ld Patient Cable,3 metri AHA/IEC, Round 11p>Din 3LD, no resistance fit for Din style LeadWires		154		101
			KV-02010X020	3Ld Patient Cable,3 metri IEC, Round 11p>NK 6J no resistance fit for BR-019 LeadWires		154		155
Nihon Kohden	compatibile BSM-2301, BSM-2353, BSM-5100, Multi-parameter bedside monitors	JC-906P	KV-0202067/4B	6Ld Patient Cable,3 metri IEC, Rectangle 12p>NK 8J, no resistance fit for BR-903P, BR-906P		173		174
Nihon Kohden	compatibile OPV-1500 Lifescope N	L060 JC-103T	KV-02010X021	3Ld Patient Cable,3 metri IEC, Rectangle 12p>NK 8J, no resistance fit for BR-903P, BR-906P				174
PHILIPS	compatibile 78352C,7834C, 78834C, M1001A/B, M1002A/B, M1722A/B, M1723A/B, M2475B, M3000A, M3001A, M4735A, VM4/V24	M1500A	KV-02010X022	3Ld Patient Cable,3 metri AHA, Round 12p>AA-3LD 1KQ resistance fit for AA style LeadWires		018		086
		M1510A	KV-02010029	3Ld Patient Cable,3 metri IEC, Round 12p>AA-3LD 1KQ resistance fit for AA style LeadWires		018		086
		M1520A	KV-02020X021	5Ld Patient Cable,3 metri AHA, Round 12p>AA-5LD 1KQ resistance fit for AA style LeadWires		018		085
		M1530A	KV-02020045	5Ld Patient Cable,3 metri IEC, Round 12p>AA-5LD 1KQ resistance fit for AA style LeadWires		018		085
PHILIPS	compatibile M3000A, M3001A, M1001A/B, M1002A/B, 78352C, 78354C	M1669A	KV-0201059/E	3Ld Patient Cable,3 metri IEC/AAMI, Round 12p>3 vie 1KQ resistance fit for Philips LeadWires		519		
		M1668A	KV-0202059/F	5Ld Patient Cable,3 metri IEC/AAMI, Round 12p>5 vie 1KQ resistance fit for Philips LeadWires		519		
	compatibile Philips Instruments connettore 8-Pin: 78208A,78210A,78330A,78331A, 78332A,78333A,78341A,78342A, 78345A,78346A,78351T,78352A,7 8353A,78354A,78532A/B,78534A/B/C, M1722A/B,M3500B	M1600A	KV-02020X022	5Ld Patient Cable,3 metri AHA, Round 8p>AA-5LD 1KQ resistance fit for AA style LeadWires		019		085
		M1610A	KV-0202034	5Ld Patient Cable,3 metri IEC, Round 8p>AA-5LD 1KQ resistance fit for AA style LeadWires		019		085

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI

PHILIPS	compatibile Philips Instruments connettore 8-Pin: 77000AC,77020A/AC,77020AR,77025A, 77030A,77035A,77040A,78208A,78210A, 78330A,78331A,78332A,78333A,78341A, 78342A,78345A,78346A,78351T,78352A, 78353A/B,78354A,78532A/B,78534A/B/C, M1722A/B,M1723A/B,M2400A,M2406A,M2475B,	M1580A	KV-02010X023	3Ld Patient Cable,3 metri AHA, Round 8p>AA-3LD 1KQ resistance fit for AA style LeadWires			086
		M1590	KV-02010024	3Ld Patient Cable,3 metri IEC, Round 8p>AA-3LD 1KQ resistance fit for AA style LeadWires			086
		M1620A	KV-02010X024	3Ld Patient Cable,3 metri AHA, Round 5p>AA-3LD 1KQ resistance fit for AA style LeadWires			086
			KV-0201027	3Ld Patient Cable,3 metri IEC, Round 5p>AA-3LD 1KQ resistance fit for AA style LeadWires			086
S&W (Artema)	compatibile 8000 System, Diascope1,2 , Neoscope, DMS 700, defibrillator, Diascope-NT, Athena Systems, Cardio-Aid, DS Traveller, DSNT, DS1,DS2, MCX/LS		KV-02020X023	5Ld Patient Cable,3 metri AHA/IEC, Round 10I>LL 5LD, no resistance fit for LL style LeadWires			105
SPACELABS	compatibile 90496 Ultraview, 90369 Ultraview, 90367 Ultraview	IEC:700-0008-07AHA:700-0008-06	KV-02020016	5Ld Patient Cable,3 metri AHA/IEC, Round 17p>TRU-LINK 5Ld, 1KQ resistance fit for TRU-LINK LeadWires			091
Fukuda Denshi	compatibile FUKUDA: 7100 patient monitor		KV-02010X025	3Ld Patient Cable,3 metri AHA/IEC, Round 12P>Euro 3LD, fit for Euro style LeadWires			242
			KV-02020X024	5Ld Patient Cable,3 metri AHA/IEC, Round 12P>Euro 5LD, fit for Euro style LeadWires			243
	compatibile DATASCOPE:Expert FUKUDA DENSHI:DS5100E,DS5300W,DS5400, DS7101/L,DS7141		KV-0202040/1	5Ld and 3Ld Patient Cable,3 metri AHA/IEC, Round 12pin> 16 jack non- resistance fit for EY032531 LeadWires			373
Lohmeier			KV-02010X026	3Ld Patient Cable,3 metri IEC/AAMI , Spin-5Jack plug>Din-3LD non-resistance fit for Din1.5 style Leads			101
	compatibile D-81241, original P/N.: 98.060.120/A	98.060.120/A	KV-02020X025	5Ld Multi-Cable,3 metri IEC/AAMI , 27pin plug>Euro 5Ld ECG+2*IBP+SpO2+2*Temp (discharge inside)			







LEAD WIRES (TERMINALI)

Compatible Brand	Productor Guide	OEM P/N	TH CODE	DESCRIPTION (Name,Length,Leads)	Color coding / Leads	Patient end termination
Bionet	compatibile con tutti gli LL-Plug System , Yoke type: Lifelinc(LL), Code: B		KV-S02010X009	LL lead wires, 24 in. (61 cm) , 3 vie	AHA	SNAP (bottone)
Datex			KV-S02010X010		IEC	CLIP (pinzetta)
Siemens			KV-S02020X009	LL lead wires, 24 in. (61 cm) , 5 vie	AHA	SNAP (bottone)
GE-Hellige			KV-S02020X010		AHA	CLIP (pinzetta)

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
il presente listino può essere soggetto a variazioni senza preavviso





LISTINO PREZZI RICAMBI COMPATIBILI

Kontron Goldway S&W (Artem)			KV-S02010X011	LL lead wires,	AHA	SNAP (bottone)	
			KV-S02010X012	40 in. (102 cm) , 3 vie	AHA	CLIP (pinzetta)	
			KV-S02010X013		IEC	SNAP (bottone)	
			KV-S02010X014		IEC	CLIP (pinzetta)	
			KV-S02020X011	LL lead wires,	AHA	SNAP (bottone)	
			KV-S02020X012	40 in. (102 cm) , 5 vie	AHA	CLIP (pinzetta)	
			KV-S02020X013		IEC	SNAP (bottone)	
			KV-S02020X014		IEC	CLIP (pinzetta)	
			KV-S02010X015	Din lead wires,	AHA	SNAP (bottone)	
			KV-S02010X016	18 in. (45 cm) , 3 vie	IEC	SNAP (bottone)	
			KV-S02020X015	Din lead wires,	AHA	SNAP (bottone)	
			KV-S02020X016	18 in. (45 cm) , 5 vie	AHA	CLIP (pinzetta)	
			KV-S02020X017		IEC	CLIP (pinzetta)	
			KV-S02010X017	Din lead wires,	AHA	SNAP (bottone)	
	KV-S02010X018	24 in. (61 cm) , 3 vie	IEC	CLIP (pinzetta)			
	KV-S02020X018	Din lead wires,	AHA	SNAP (bottone)			
	KV-S02020X019	24 in. (61 cm) , 5 vie	AHA	CLIP (pinzetta)			
	KV-S02020X020		IEC	SNAP (bottone)			
	KV-S02020X021		IEC	CLIP (pinzetta)			
	KV-S02010X019	Din lead wires,	AHA	SNAP (bottone)			
	KV-S02010X020	40 in. (102 cm) , 3 vie	AHA	CLIP (pinzetta)			
	KV-S02010X021		IEC	SNAP (bottone)			
	KV-S02010X022		IEC	CLIP (pinzetta)			
	KV-S02020X022	Din lead wires,	AHA	SNAP (bottone)			
	KV-S02020X023	40 in. (102 cm) , 5 vie	AHA	CLIP (pinzetta)			
	KV-S02020X024		IEC	SNAP (bottone)			
	KV-S02020X025		IEC	CLIP (pinzetta)			
Biosys Creative Datascope Nihon Kohden AIR SHIELDS (NARCO, Vickers) ATL ULTRASOUND	compatibile con yoke type: DIN:	  					
			M1673A	KV-S02010X023	Multi-Link Ldwr,3 vie	AHA	SNAP (bottone)
			M1674A	KV-S0201059/A	Set-Group, 40 in (102cm)	IEC	SNAP (bottone)
			M1671A	KV-S02010X024		AHA	CLIP (pinzetta)
			M1672A	KV-S0201059/B		IEC	CLIP (pinzetta)
			M1644A,	KV-S02020X026	Multi-Link Ldwr,5 vie	AHA	SNAP (bottone)
			M1645A,	KV-S0202059/C	Set-Group,	IEC	SNAP (bottone)
			M1968A,	KV-S02020X027	LA,RA,V=40in(102cm);	AHA	CLIP (pinzetta)
			M1971A,	KV-S0202059/D	LL,RL=61in(155cm);	IEC	CLIP (pinzetta)
			M1605A	KV-S02010X025	AA lead wires,	AHA	SNAP (bottone)
		M1603A	KV-S02010X026	40 in. (102 cm) , 3 vie	AHA	CLIP (pinzetta)	
PHILIPS	compatibile con tutti i AA-system Yoke type: AA style, Code: D, Or: per utilizzo in sala operatoria; include filtro specifico; non consente il monitoraggio respiratorio		M1615A	KV-S02010054		IEC	SNAP (bottone)
		M1613A	KV-S02010056		IEC	CLIP (pinzetta)	
		M1625A	KV-S02020X028	AA lead wires, 5 vie	AHA	SNAP (bottone)	
		M1623A	KV-S02020X029	LA,RA,V 40in(102cm);	AHA	CLIP (pinzetta)	
		M1635A	KV-S02020055/I	LL,RL 61in(155cm)	IEC	SNAP (bottone)	
							

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso














LISTINO PREZZI RICAMBI COMPATIBILI

		M1633A	KV-S0202059		IEC	CLIP (pinzetta)
		M1601A	KV-S02010X027	AA lead wires, OR	AHA	CLIP (pinzetta)
		M1611A	KV-S0201059/C	40 in. (102 cm) , 3 vie	IEC	CLIP (pinzetta)
		M1621A	KV-S02020X030	AA lead wires, 5 vie , OR	AHA	CLIP (pinzetta)
		M1631A	KV-S0202059/B	LA,RA,V 40in(102cm); LL,RL 61in(155cm)	IEC	CLIP (pinzetta)
GE--Marquette (Set --Group type Lead wires)	compatibile con tutti i Multi-link (VS)-plug system yoke type: Multi-Link (VS), Code: E 	E9008LF	KV-S02010X028	Multi-Link Ldwr	AHA	SNAP (bottone)
		E9003CL	KV-S02010X029	Set -- Group,,	AHA	CLIP (pinzetta)
		E9003CN	KV-S020109/a	29 in. (76 cm) , 3 vie	IEC	CLIP (pinzetta)
		E9008LG	KV-S02010X030	Multi-Link Ldwr	AHA	SNAP (bottone)
		E9003CM	KV-S02010X031	Set -- Group,,	AHA	CLIP (pinzetta)
		E9003CP	KV-S020109/a1	51 in. (130 cm) , 3 vie	IEC	CLIP (pinzetta)
GE--Marquette (Set --Group type Lead wires)	compatibile con tutti i Multi-link (VS)-plug system yoke type: Multi-Link (VS), Code: E 	E9008L	KV-S02020X031	Multi-Link Ldwr	AHA	SNAP (bottone)
		E9003CE	KV-S02020X032	Set -- Group,,	AHA	CLIP (pinzetta)
		E9008LB	KV-S0203012	29 in. (76 cm) , 5 vie	IEC	SNAP (bottone)
		E9003CG	KV-S02020X033		IEC	CLIP (pinzetta)
		E9003CF	KV-S02020X034	Multi-Link Ldwr	AHA	CLIP (pinzetta)
			KV-S02020X035	Set -- Group,,	IEC	CLIP (pinzetta)
		E9008LD	KV-S02020X036	Multi-Link Leadwire Set Group, 5 vie , Mixed Length, R _L ,L _L ,V--51in.(130cm) ; R _A ,L _A --29 in.(74cm)	AHA	SNAP (bottone)
		E9003CJ	KV-S02020X037		AHA	CLIP (pinzetta)
			KV-S02020X038	Multi-Link Leadwire Set Group, 5 vie , Mixed Length, N,F,C--51in. (130cm) , R _L --29 in.(74cm)	IEC	CLIP (pinzetta)
GE-Medical Marquette (Set --Assembly type Lead wires)	compatibile con tutti i Multi-link (VS)-plug system yoke type: Multi-Link (VS), Code: E 	E9008KB	KV-S02020X039	Multi-Link Leadwire Set	AHA	SNAP (bottone)
		E9002ZW	KV-S02020X040	29 in. (76 cm) , 5 vie	AHA	CLIP (pinzetta)
		E9008KD	KV-S0202012		IEC	SNAP (bottone)
		E9002ZZ	KV-S0203007		IEC	CLIP (pinzetta)
		E9008KC	KV-S02020X041	Multi-Link Leadwire Set ,	AHA	SNAP (bottone)
		E9002ZY	KV-S02020X042	51 in. (130 cm) , 5 vie	AHA	CLIP (pinzetta)
		E9008KE	KV-S02020X043		IEC	SNAP (bottone)
		E9003A	KV-S020208		IEC	CLIP (pinzetta)
		E9008KF	KV-S02020X044	Multi-Link Leadwire Set , 5 vie , Mixed Length, R _L ,L _L ,V--51in.(130cm); R _A ,L _A --29 in.(74cm)	AHA	SNAP (bottone)
		E9003AA	KV-S02020X045		AHA	CLIP (pinzetta)
			KV-S02020X046	Multi-Link Leadwire Set , 5 vie , Mixed Length, N,F,C--51in.(130cm) ; R _L --29 in.(74cm)	IEC	SNAP (bottone)
		E9003AB	KV-S020209		IEC	CLIP (pinzetta)
Datex PHILIPS PIONEER PPG SPACELABS BURDICK TEKTRONIX (Vitatek) WELCH ALLYN	compatibile con tutti i CONMED-plug system yoke type: CONMED, Code: G 		KV-S02010X032	Conmed lead wires, 24 in. (61 cm) , 3 vie	IEC	CLIP (pinzetta)
			KV-S02020X047	Conmed lead wires,	AHA	SNAP (bottone)
			KV-S02020X048	24 in. (61 cm) , 5 vie	AHA	CLIP (pinzetta)
			KV-S02020X049		IEC	CLIP (pinzetta)
			KV-S02010X033	Conmed lead wires,	AHA	SNAP (bottone)
			KV-S02010X034	40 in. (102 cm) , 3 vie	AHA	CLIP (pinzetta)
			KV-S02010X035		IEC	SNAP (bottone)
			KV-S02010X036		IEC	CLIP (pinzetta)
			KV-S02020X050	Conmed lead wires,	AHA	SNAP (bottone)
			KV-S02020X051	40 in. (102 cm) , 5 vie	AHA	CLIP (pinzetta)
			KV-S02020X052		IEC	SNAP (bottone)
			KV-S02020X053		IEC	CLIP (pinzetta)

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
il presente listino può essere soggetto a variazioni senza preavviso





LISTINO PREZZI RICAMBI COMPATIBILI

Datex-Ohmeda GE-medical Pro1000 	compatibile con tutti gli ASP-plug system yoke type: ASP, Code: H  EH024C3I	KV-S02010X037	ASP lead wires, 24 in. (61 cm) , 3 vie	AHA	SNAP (bottone)		
		KV-S02010X038	IEC	CLIP (pinzetta)			
		KV-S02020X054	ASP lead wires, 24 in. (61 cm) , 5 vie	IEC	CLIP (pinzetta)		
		KV-S02010X039	ASP lead wires, 40 in. (102 cm) , 3 vie	AHA	SNAP (bottone)		
		KV-S02010X040		AHA	CLIP (pinzetta)		
		KV-S0202064/1a		IEC	SNAP (bottone)		
		KV-S02010064/1A	ASP lead wires, 40 in. (102 cm) , 5 vie	IEC	CLIP (pinzetta)		
		KV-S02020X055		AHA	SNAP (bottone)		
		KV-S02020X056		AHA	CLIP (pinzetta)		
		KV-S0202064/1		IEC	SNAP (bottone)		
KV-S02010064/1	IEC	CLIP (pinzetta)					
Type LM Leadwires 494 	 	KV-S02010X042	LM lead wires, 40 in. (102 cm) , 3 vie	AHA	SNAP (bottone)		
		KV-S02010X043	AHA	CLIP (pinzetta)			
		KV-S02010X044	IEC	SNAP (bottone)			
		KV-S02010X045	IEC	CLIP (pinzetta)			
		KV-S02020X058	LM lead wires, 40 in. (102 cm) , 5 vie	AHA	SNAP (bottone)		
		KV-S02020X059	AHA	CLIP (pinzetta)			
		KV-S02020X060	IEC	SNAP (bottone)			
		KV-S02020X061	IEC	CLIP (pinzetta)			
		Mindray Siemens SPACELABS BURDICK Viggo- Spectramed (Statham) WESTINGHOUSE 084 	compatibile con tutti gli EURO-plug system yoke type: EURO-style, Code: J  EJ040S3I  EJ040S5I	KV-S02010X046	Euro style lead wires, 24 in. (61 cm) , 3 vie	IEC	SNAP (bottone)
				KV-S02010X047	IEC	CLIP (pinzetta)	
KV-S02020X062	Euro style lead wires, 24 in. (61 cm) , 5 vie			AHA	CLIP (pinzetta)		
KV-S02020X063				IEC	CLIP (pinzetta)		
KV-S02010X048				AHA	SNAP (bottone)		
KV-S02010X049	Euro style lead wires, 40 in. (102 cm) , 3 vie			AHA	CLIP (pinzetta)		
KV-S02010X050	AHA			Mini Grabber (Clip)			
KV-S0102021	IEC			SNAP (bottone)			
KV-S010201	IEC			CLIP (pinzetta)			
KV-S02020X064	Euro style lead wires, 40 in. (102 cm) , 5 vie			AHA	SNAP (bottone)		
KV-S02020X065	AHA	CLIP (pinzetta)					
KV-S0202091/8A	IEC	SNAP (bottone)					
KV-S020301	IEC	CLIP (pinzetta)					
NEC suitable for all NEC-Plug system, Yoke type: NEC style code: K  097  100	47505	KV-S02010X051	47504 style lead wires, 29 in. (76 cm) , 3 vie	IEC	SNAP (bottone)		
		KV-S02020X067	47505 style lead wires, 29 in. (76 cm) , 5 vie	IEC	SNAP (bottone)		
		KV-S02020X068	IEC	CLIP (pinzetta)			
SPACELABS 166 	compatibile con tutti i TRU-LINK-plug system yoke type: TRU-LINK style, Code: L  EL040S3I  EL040S5I	700-0007-00	TUR-LINK style lead wires, 24 in. (61 cm) , 3 vie	AHA	SNAP (bottone)		
		700-0007-04		KV-S02010X052/A	AHA	SNAP (bottone)	
		700-0006-04		KV-S02010X052/B	AHA	CLIP (pinzetta)	
		700-0007-05	KV-S02010X052/C	IEC	SNAP (bottone)		
		700-0006-05	KV-S02010X052/D	IEC	CLIP (pinzetta)		
		700-0007-08	KV-S02020X068	TUR-LINK style lead wires, 24 in. (61 cm) , 5 vie	AHA	SNAP (bottone)	
		700-0007-09	KV-S02020X069		IEC	SNAP (bottone)	
		700-0006-09	KV-S02020X070		IEC	CLIP (pinzetta)	
		700-0007-12	KV-S02020X070/A	TUR-LINK style lead wires, 24 in. (61 cm) , 6 vie	AHA	SNAP (bottone)	
		700-0006-12	KV-S02020X070/B		AHA	CLIP (pinzetta)	
		700-0007-13	KV-S02020X070/C		IEC	SNAP (bottone)	
		700-0006-13	KV-S02020X070/D	IEC	CLIP (pinzetta)		
		700-0007-16	KV-S03020X001	TUR-LINK style lead wires, 24 in. (61 cm) , 10 vie	AHA	SNAP (bottone)	
		700-0006-16	KV-S03020X002		AHA	CLIP (pinzetta)	
		700-0007-17	KV-S03020X003		IEC	SNAP (bottone)	
		700-0006-17	KV-S03020X004		IEC	CLIP (pinzetta)	

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso





LISTINO PREZZI RICAMBI COMPATIBILI

		700-0007-03	KV-S02010X053	TUR-LINK style lead wires, 40 in. (102 cm) , 3 vie	IEC	SNAP (bottone)
		700-0006-03	KV-S0201020		IEC	CLIP (pinzetta)
		700-0007-06	KV-S02010X053/A		AHA	SNAP (bottone)
		700-0006-06	KV-S02010X053/B		AHA	CLIP (pinzetta)
		700-0007-07	KV-S02010X053/C	TUR-LINK style lead wires, 40 in. (102 cm) , 4 vie	IEC	SNAP (bottone)
		700-0006-07	KV-S02010X053/D		IEC	CLIP (pinzetta)
		700-0007-10	KV-S02020X071		AHA	SNAP (bottone)
		700-0006-10	KV-S02020X072	TUR-LINK style lead wires, 40 in. (102 cm) , 5 vie	AHA	CLIP (pinzetta)
		700-0007-11	KV-S0202023		IEC	SNAP (bottone)
		700-0006-11	KV-S0202021		IEC	CLIP (pinzetta)
		700-0007-14	KV-S02020X072/A		AHA	SNAP (bottone)
		700-0006-14	KV-S02020X072/B		AHA	CLIP (pinzetta)
		700-0007-15	KV-S02020X072/C	TUR-LINK style lead wires, 40 in. (102 cm) , 6 vie	IEC	SNAP (bottone)
		700-0006-15	KV-S02020X072/D		IEC	CLIP (pinzetta)
		700-0007-18	KV-S03020X005		AHA	SNAP (bottone)
		700-0006-18	KV-S03020X006		AHA	CLIP (pinzetta)
		700-0007-19	KV-S03020X007	TUR-LINK style lead wires, 40in. (102 cm) , 10 vie	IEC	SNAP (bottone)
		700-0006-19	KV-S03020X008		IEC	CLIP (pinzetta)
Colin compatibile con tutti i Colin- Plug System , Yoke type: Colin style, Code: M	 009	CL-3601	KV-S02010X054	COLIN style lead wires, 40in. (102 cm) , 3 vie	IEC	SNAP (bottone)
			KV-S02010X055		IEC	CLIP (pinzetta)
	 010		KV-S02020X073		COLIN style lead wires, 40in. (102 cm) , 5 vie	IEC
Nihon Konden compatibile con tutti i NK- 5LD Plug System , Yoke type: NK-5LD style, Code: N	 126	BR-019	KV-S0201067/8	NK style lead wires, 32 in. 80 cm , 3 vie	IEC	SNAP (bottone)
			KV-S0201067/3		IEC	CLIP (pinzetta)
		BR-019P	KV-S02020X074		IEC	SNAP (bottone)
	 125	BR-021P	KV-S0202067/2		NK style lead wires, 32 in. 80 cm , 5 vie	IEC





LISTINO CAVI ECG COMPATIBILI



PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso



LISTINO PREZZI RICAMBI COMPATIBILI

Nihon Konden compatibile con tutti i NK-6LD Plug System , Yoke type: NK-6LD style, Code: P	 171/433	BR-913P	KV-S02010X056	NK style lead wires, 35 in. 90 cm , 3 vie	IEC	SNAP (bottone)
		BR-903P	KV-S0201067/1A		IEC	CLIP (pinzetta)
			KV-S02020X074/A	NK style lead wires, 35 in. 90 cm , 6 vie	IEC	SNAP (bottone)
		BR-906P	KV-S0202067/1		IEC	CLIP (pinzetta)
Nihon Konden compatibile con tutti il NK-5LD Plug System , Yoke type: NK-5LD style, Code: R	 124		KV-S0201067/7A	NK style lead wires, 32 in. 80 cm , 3 vie	IEC	SNAP (bottone)
		BR-546S	KV-S0201067/7		IEC	CLIP (pinzetta)
	 123		KV-S0201067/5	NK style lead wires, 32 in. 80 cm , 3 vie	IEC	CLIP (pinzetta)
		BR-002P	KV-S02020X074/A		IEC	SNAP (bottone)
FUKUDA DENSHI  372		KV-S02010X057	one-piece lead wires,32 in. 80 cm , 3 vie	IEC	SNAP (bottone)	
		KV-S02010X058		AHA	SNAP (bottone)	

Compatible Philips New ECG Monitoring Cable and Lead System


Compatible Brand	Compatible Models N.	OEM P/N	TH CODE	DESCRIPTION (Name,Length,Plug,Compatible)	instrument connector	YOKE (lead side)
PHILIPS	compatibile M3000ACM3001AC M1001A/B,M1002A/B, 78352C,78354C	M1669A	KV-0201059/E	cavo paziente 3 vie (cavo base) , 2,7 metri IEC/AAMI, Round 12p>3 vie 1KΩ resistance fit for Philips LeadWires		
	compatibile M3000A, M3001A, M1001A/B,M1002A/B, 78352C,78354C	M1668A	KV-0202059/F	cavo paziente 5 vie (cavo base) , 2,7 metri IEC/AAMI, Round 12p>5 vie 1KΩ resistance fit for Philips LeadWires		

Compatible Brand	Productor Guide	OEM P/N	TH CODE	DESCRIPTION (Name,Length,Leads)	Color coding / Leads	Patient end termination
PHILIPS Compatibile M1668A 5-lead trunk cable. Use EDG073-3AI or EDG073-5AI Trunk cable.		M1673A	KV-S02010X059	Multi-Link Ldwr, terminali 3 vie Set-Group, 40 in (102cm)	AHA	SNAP (bottone)
		M1674A	KV-S0201059/A		IEC	SNAP (bottone)
		M1671A	KV-S02010X060		AHA	CLIP (pinzetta)
		M1672A	KV-S0201059/B		IEC	CLIP (pinzetta)
		M1644A,	KV-S02020X075	Multi-Link Ldwr, terminali 5 vie Set-Group, LA,RA,V=40in(102cm); LL,RL=61in(155cm);	AHA	SNAP (bottone)
		M1645A,	KV-S0202059/C		IEC	SNAP (bottone)
		M1968A,	KV-S02020X076		AHA	CLIP (pinzetta)
		M1971A,	KV-S0202059/D		IEC	CLIP (pinzetta)


Compatible Brand	Compatible Models N.	OEM P/N	TH CODE	DESCRIPTION (Name,Length,Plug,Compatible)	instrument connector	YOKE (lead side)
PHILIPS	compatibile 78352C, 78354C, 78834C Legacy HP Monitors,Modules M1001A/B,M1002A/B, M3001A Opt. C12 or Opt. C18 IntelliVue Multi-Measurement Modules M3002A Opt. C12 MMS X2 M8102A Opt. C12 MP2 Monitor M8105A Opt. C12 MP5 Monitor M3535A/M3536A HeartStar	M1949A	KV-0302X186	ECG trunk cable (cavo base) connettore 12-pin. AAMI and IEC labels included. Length: 8.9' (2.7m)		




LISTINO PREZZI RICAMBI COMPATIBILI

Compatible Brand	Productor Guide	OEM P/N	TH CODE	DESCRIPTION (Name,Length,Leads)	Color coding / Leads	Patient end termination
PHILIPS Compatibile M1663A or M1949A		M1644A+ M1602A	KV-S03020X009	5+5 Lead Sets, Length:1m, Lead mark: AAMI: (RL,LA,LL,RA,V) + (V2,V3,V4,V5,V6;) IEC: (N, L, F, R, C) + (C2,C3,C4,C5,C6;)	AHA	SNAP (bottone)
		M1645A+ M1604A	KV-S03020X010		IEC	SNAP (bottone)
		M1968A+ M1976A	KV-S03020X011		AHA	CLIP (pinzetta)
		M1971A+ M1978A	KV-S03020X012		IEC	CLIP (pinzetta)

One-Piece Series EKG Cable With Leads (Cavo completo 10 vie con terminali pressofusi)










Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	TH CODE	STD	PATIENT END TERMINAL	defibrillation withstand
Bosch Cambridge Custo-med DEGO Hellige Siemens Von Berg (hormann)	compatibile: Bosch:603 D, 603 ER, 603 S, 606		KV-0302X001	AHA	CLIP	No Resistance
	Cambridge: CM 6000, CM 6400		KV-0302X002	AHA	DIN3.0	No Resistance
	Custo-med: Custo-card		KV-0302X003	AHA	BANANA	No Resistance
	DEGO: Cardiomed 3		KV-0302X004	AHA	SNAP	4.7kΩ Resistance
	Hellige: EK 36, EK 41, EK 43, EK 53,		KV-0302073	IEC	SNAP	4.7kΩ Resistance
	EK 56, EK 403, EK 413, EK 512, EK 512 E,		KV-0302X005	AHA	DIN3.0	4.7kΩ Resistance
	EK 512 P, Cardiosmart, CardioSys, Corna		KV-0302073/1	IEC	CLIP	No Resistance
	Siemens:Cardiostat1/3T/703/710/31/31S		KV-0302X006	IEC	DIN3.0	No Resistance
	Cardina C, Cardirex 3T/4/6T/62, Cathcor, Megacard Elema A/B, R, R/E, Mingodoc, Mingocard 3 / 4 / 7, Mingograph 4 / 62 / 81 Mingograph 82 / 740, Mingolog 7,		KV-0302X007	IEC	BANANA	No Resistance
	Mingophon 3/7/310/410/420/710/720, Miracard 1, Sirecard440/440S/460/460S/740		KV-0302X008	IEC	DIN3.0	4.7kΩ Resistance
	Siereg, Sirecard 803, Sirecust 701		KV-0302X009			
	Von Berg:Bioset3500/3600/3700/6000/8000		KV-0302072	IEC	BANANA	4.7kΩ Resistance

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI

Hellige	compatibile Hellige: EK 51		KV-0302X010	AHA	CLIP	10 kΩ Resistance	
			KV-0302X011	IEC	DIN3.0	10 kΩ Resistance	
Avionics(Del Mar) Custo-med Del mar Ergoline Esaote Biomed Fidelity Innomed Kontron Reynolds Schiller Trentina (Elettronica) Von Berg Bionet Welch allyn,	compatibile: Avionics(Del Mar): 910 / 920 / 930 / 940 / 942 / 950 / 960 / 970 / 980 Custo-med: Custo-norm Del mar: 910,920,930,940,942,950,960,970,980 Ergoline: Autoruler, Autoscript 6/12 Cardiette, EK 3003 / 3012, Ergoline Esaote Biomedica: EKG P120 Fidelity: LP 3000 Innomed: Heart Mirror (3-Kanal) Kontron: EKG P80 Reynolds: Cardiotrack 12 Schiller:AT3, AT5, AT6,AT10, AT60,AT104 PC, AT104, CS6, CS12,CS100, CS200 Trentina: All Excel-Types Von Berg: Bioset 9000		KV-0302X012	AHA	SNAP	No Resistance	
			KV-0302X013	AHA	BANANA	No Resistance	
			KV-0302X014	AHA	BANANA	10 kΩ Resistance	
			KV-0302X015	AHA	BANANA	10 kΩ+discharge tube	
			KV-0302X016	IEC	CLIP	No Resistance	
			KV-0302X017	IEC	BANANA	No Resistance	
			KV-03020176	IEC	CLIP	10 kΩ Resistance	
			KV-0302X018	IEC	DIN3.0	10 kΩ Resistance	
			KV-0302080	IEC	BANANA	10 kΩ Resistance	
			Schiller Esaote Biomedica	compatibile Schiller: AT-2 Plus, AT1, AT2, AT4, AT101, AT102 Esaote Biomedica: P80, P80Power, P80PowerSix Mortara ELI 250 ELI 350		KV-0302X019	AHA
KV-0302X020	AHA	CLIP				10 kΩ Resistance	
KV-0302X021	AHA	DIN3.0				10 kΩ Resistance	
KV-0302X022	AHA	BANANA				10 kΩ Resistance	
KV-0302083	IEC	SNAP				10 kΩ Resistance	
KV-0302083/1	IEC	CLIP				10 kΩ Resistance	
KV-0302X025	IEC	DIN3.0				10 kΩ Resistance	
KV-0302082	IEC	BANANA				10 kΩ Resistance	
Bosch / dimeq Cambridge Cardioline/Remco Cardionics Hellige Honeywell/Philips: Polymed Schwarzer SECA	compatibile Bosch/dimeq: 501A,501NV,503,601 503B,503D,503NV,503V,506,506D, Cambridge: VS 6, SC 1000 Delta 1 / 3 / 300 / 360, E 1 / E 1 R / E 3 Cardionics: all Models Hellige: Corina. EK 51 Honeywell/Philips: Cardioplan, CR 300, CR 330, ER 330 Polymed: Alpha 600 / 800 Schwarzer: VS 6 SECA: 3000		KV-0302X027	AHA	SNAP	4.7kΩ Resistance	
			KV-0302X028	AHA	CLIP	4.7kΩ Resistance	
			KV-0302X029	AHA	DIN3.0	4.7kΩ Resistance	
			KV-0302X030	AHA	BANANA	4.7kΩ Resistance	
			KV-0302X031	IEC	SNAP	4.7kΩ Resistance	
			KV-0302X032	IEC	CLIP	4.7kΩ Resistance	
			KV-0302X033	IEC	DIN3.0	4.7kΩ Resistance	
			KV-0302073/4	IEC	BANANA	4.7kΩ Resistance	
Burdick Quinton ATRIA	compatibile BURDICK (QUINTON):Eclipse 4, Eclipse 400,Eclipse 8, Eclipse 800, Eclipse 850, Eclipse 850i, Eclipse Plus, Eclipse LE II BURDICK with serial number starting with 1067 or 92304 ATRIA: ATRIA PLUS, ATRIA3000		KV-0302X035	AHA	CLIP	10 kΩ Resistance	
			KV-0302X036	AHA	DIN3.0	10 kΩ Resistance	
			KV-0302X037	AHA	BANANA	10 kΩ Resistance	
			KV-0302X038	IEC	CLIP	10 kΩ Resistance	
			KV-0302X039	IEC	DIN3.0	10 kΩ Resistance	
			KV-0301073/5	IEC	BANANA	10 kΩ Resistance	
Burdick	compatibile EK-10		KV-0302X040	AHA	CLIP	No Resistance	
			KV-0302X041	AHA	BANANA	No Resistance	
			KV-0302X042	IEC	SNAP	No Resistance	
			KV-0302X043	IEC	CLIP	No Resistance	
	compatibile Burdick: Elite, Elite II, E-350, E350i, E550, E560, (1st Generation Eclipse LE/LEi, with serial number starting with 1051).		KV-0302X044	IEC	DIN3.0	No Resistance	
			KV-0302X045	IEC	BANANA	No Resistance	
			KV-0302X046	AHA	SNAP	4.7kΩ Resistance	
			KV-0302X047	AHA	CLIP	4.7kΩ Resistance	

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI

	Stress cable with detachable leads 007706 Fukuda Denshi: Alpha 500 / 600, FCP 2201		KV-0302X048	AHA	DIN3.0	4.7kΩ Resistance		
			KV-0302X049	AHA	BANANA	4.7kΩ Resistance		
			KV-0302X050	IEC	SNAP	4.7kΩ Resistance		
EDAN	compatibile SE-12 Express SE-3 SE-601A		KV-0302X051	AHA	SNAP	10 kΩ Resistance		
			KV-0302X052	AHA	CLIP	10 kΩ Resistance		
			KV-0302X053	AHA	DIN3.0	10 kΩ Resistance		
			KV-0302X054	AHA	BANANA	10 kΩ Resistance		
			KV-0302X055	IEC	CLIP	10 kΩ Resistance		
			KV-0302X056	IEC	DIN3.0	10 kΩ Resistance		
			KV-0302X057	IEC	BANANA	10 kΩ Resistance		
						KV-0302X058	AHA	BANANA
Fukuda Denshi	compatibile Fukuda Denshi Autocardiner, Cardimax FX-2111 /FX-3010 / FX-4010, FCP-2155		KV-0302X059	IEC	DIN3.0	4.7kΩ Resistance		
			KV-0302X060	IEC	BANANA	4.7kΩ Resistance		
			KV-0302X061	AHA	CLIP	No Resistance		
FUKUDA ME	compatibile KP-500		KV-0302X062	AHA	BANANA	No Resistance		
			KV-0302X063	IEC	SNAP	No Resistance		
			KV-0302X064	IEC	CLIP	No Resistance		
			KV-0302X065	IEC	DIN3.0	No Resistance		
	compatibile KP-500D			KV-0302X066	AHA	SNAP	20kΩ+discharge	
				KV-0302X067	AHA	CLIP	20kΩ+discharge	
				KV-0302X068	AHA	DIN3.0	20kΩ+discharge	
				KV-0302X069	IEC	SNAP	20kΩ+discharge	
				KV-0302X070	IEC	CLIP	20kΩ+discharge	
				KV-0302X071	IEC	DIN3.0	20kΩ+discharge	
	compatibile FX-031			KV-0302X072	AHA	SNAP	No Resistance	
				KV-0302X073	AHA	CLIP	No Resistance	
				KV-0302X074	IEC	SNAP	No Resistance	
				KV-0302X075	IEC	CLIP	No Resistance	
				KV-0302X076	IEC	DIN3.0	No Resistance	
				KV-0302X077	AHA	SNAP	4.7kΩ Resistance	
compatibile Fukuda M-E: Cardisuny alpha 1000, Cardisuny alpha 6600 A / 6600 AX			KV-0302X078	AHA	CLIP	4.7kΩ Resistance		
			KV-0302X079	AHA	DIN3.0	4.7kΩ Resistance		
			KV-0302X080	AHA	BANANA	4.7kΩ Resistance		
			KV-0302X081	IEC	SNAP	4.7kΩ Resistance		
			KV-0302X082	IEC	CLIP	4.7kΩ Resistance		
			KV-0302X083	IEC	DIN3.0	4.7kΩ Resistance		
			KV-0302X084	IEC	BANANA	4.7kΩ Resistance		
			KV-0302X085	AHA	CLIP	No Resistance		
			KV-0302X086	AHA	DIN3.0	No Resistance		
			KV-0302X087	AHA	BANANA	No Resistance		
KANZ PC-104, Cardioline remco	Kanz: 103, 106 Cardioline/remco: delta 1,3,30,60, delta 1 plus, delta 3 plus,delta 60 plus		KV-0302X088	AHA	SNAP	4.7kΩ Resistance		
			KV-0302X089	AHA	CLIP	4.7kΩ Resistance		
			KV-0302X090	IEC	SNAP	No Resistance		
			KV-0302X091	IEC	DIN3.0	No Resistance		
			KV-0302X092	IEC	BANANA	No Resistance		
			KV-03020301	IEC	CLIP	4.7kΩ Resistance		
			KV-03020201	IEC	CLIP	4.7kΩ Resistance		
			KV-03020101	IEC	BANANA	4.7kΩ Resistance		

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI

Kanz PC-109	compatibile ECG 108/110/1203,1205		KV-0302X093	AHA	SNAP	No Resistance
			KV-0302X094	AHA	CLIP	No Resistance
			KV-0302X095	AHA	BANANA	No Resistance
			KV-0302X096	IEC	CLIP	No Resistance
			KV-0302X097	IEC	DIN3.0	No Resistance
			KV-0302X098	IEC	BANANA	No Resistance
MARQUETTE (GEMS)	compatibile Hellige: MicroSmart, MAC 500/1100, MAC 1200 / 1200 ST		KV-0302X099	AHA	CLIP	10 kΩ Resistance
			KV-0302X100	AHA	BANANA	10 kΩ Resistance
			KV-0302089/4	IEC	SNAP	10 kΩ Resistance
			KV-0302MQ02	IEC	BANANA	10 kΩ Resistance
			KV-0302X101	AHA	SNAP	10 kΩ Resistance
			KV-0302X102	AHA	CLIP	10 kΩ Resistance
MEDICAL RESEARCH LABS (MRL)	compatibile Portable Intensive Care (PIC) System		KV-0302X103	AHA	DIN3.0	10 kΩ Resistance
			KV-0302X104	AHA	BANANA	10 kΩ Resistance
			KV-0302069/2	IEC	SNAP	10 kΩ Resistance
			KV-0303069/2B	IEC	CLIP	10 kΩ Resistance
			KV-0302X105	IEC	DIN3.0	10 kΩ Resistance
			KV-0302X106	IEC	BANANA	10 kΩ Resistance
MEDTRONIC PHYSIO CONTROL	compatibile Lifepak 12		KV-0302X107	AHA	SNAP	10 kΩ Resistance
			KV-0302X108	AHA	CLIP	10 kΩ Resistance
			KV-0302X109	AHA	DIN3.0	10 kΩ Resistance
			KV-0301068	AHA	BANANA	10 kΩ Resistance
			KV-0301069	IEC	SNAP	10 kΩ Resistance
			KV-0303069/1	IEC	SNAP	10 kΩ Resistance
Mindray	compatibile BeneView (12-Ld module) Oem p/n:0010-30-42721 (AAMI std.) Oem p/n:0010-30-42722 (IEC std.)		KV-0302X110	IEC	DIN3.0	10 kΩ Resistance
			KV-0302X111	AHA	SNAP	10 kΩ Resistance
			KV-0302X112	AHA	CLIP	10 kΩ Resistance
			KV-0302X113	AHA	DIN3.0	10 kΩ Resistance
			KV-0302X114	AHA	BANANA	10 kΩ Resistance
			KV-0302071/1	IEC	SNAP	10 kΩ Resistance
MORTARA	compatibile MORTARA INSTRUMENT:ELI 100 MORTARA INSTRUMENT:ELI 50 MORTARA INSTRUMENT:ELI 200		KV-0303071/2	IEC	CLIP	10 kΩ Resistance
			KV-0302X115	IEC	DIN3.0	10 kΩ Resistance
			KV-0302X116	AHA	SNAP	4.7kΩ Resistance
			KV-0302X117	AHA	BANANA	4.7kΩ Resistance
			KV-0303MI	IEC	SNAP	4.7kΩ Resistance
			KV-0302MI	IEC	CLIP	4.7kΩ Resistance
MORTARA	compatibile Surveyor (New 12 pin round connector)		KV-0302X118	IEC	DIN3.0	4.7kΩ Resistance
			KV-0301MI	IEC	BANANA	4.7kΩ Resistance
NEC	compatibile NEC YY-1063		KV-0302069/3	IEC	SNAP	10 kΩ Resistance
			KV-0302X119	AHA	SNAP	10 kΩ Resistance
			KV-0302X120	AHA	CLIP	10 kΩ Resistance
			KV-0302X121	AHA	DIN3.0	10 kΩ Resistance
			KV-0302X122	AHA	BANANA	10 kΩ Resistance
			KV-0302X123	IEC	SNAP	10 kΩ Resistance
Nihon Kohden	compatibile : Nihon Konden: ECG-9620/9010/9020 ECG-9022/9110/9130/9132 ECG-8420(tranne 9320) BIOCARE邦健: ECG-101/101G/300 ECG-300G, ECG-1200, ECG-9803.		KV-0302X124	IEC	CLIP	10 kΩ Resistance
			KV-0302X125	IEC	DIN3.0	10 kΩ Resistance
			KV-0302X126	IEC	BANANA	10 kΩ Resistance
			KV-0302X127	AHA	SNAP	No Resistance
			KV-0302X128	AHA	CLIP	No Resistance
			KV-0302X129	AHA	DIN3.0	No Resistance
Biocare邦健	CMICS: ECG-11A, ECG-11B, ECG-11D, ECG-32A, ECG-32A(II), ECG-120 ,		KV-0302X130	AHA	BANANA	No Resistance
			KV-0302X131	IEC	SNAP	No Resistance
DongJiang东江	CMICS: ECG-11A, ECG-11B, ECG-11D, ECG-32A, ECG-32A(II), ECG-120 ,		KV-0302X132	IEC	CLIP	No Resistance
			KV-0302X133	IEC	DIN3.0	No Resistance
			KV-0302X134	IEC	BANANA	No Resistance

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso




LISTINO PREZZI RICAMBI COMPATIBILI

Nihon Kohden	compatibile: BSM-2301,BSM-2353, BSM-5100, Supply Code No. K901 Article No.BJ-900P		KV-0302X135	IEC	SNAP	No Resistance	
			KV-0302X136	IEC	CLIP	No Resistance	
	compatibile: Nihon Kohden: Cardiofax6151/6353, Cardiofax EKG-8270/8350/8370/8420 8423/8450/8453/8820/8823/8830/9130 Cardiofax-Q		KV-0302X137	AHA	CLIP	4.7kΩ Resistance	
			KV-0303087/1	IEC	CLIP	4.7kΩ Resistance	
			KV-0302X138	AHA	DIN3.0	4.7kΩ Resistance	
			KV-0302X139	AHA	BANANA	4.7kΩ Resistance	
			KV-03020186	IEC	SNAP	4.7kΩ Resistance	
			KV-0302X141	IEC	DIN3.0	4.7kΩ Resistance	
			KV-030289/4D	IEC	BANANA	4.7kΩ Resistance	
			KV-0302X143	AHA	SNAP	No Resistance	
			KV-0302X144	AHA	BANANA	No Resistance	
			KV-0302X145	IEC	DIN3.0	No Resistance	
	compatibile: Cardiofax Q ECG-9110K Cardiofax Q ECG-9130K, ECG-9130P, ECG-9620P special screw, imperial standard(SAE)		KV-030289/4I	IEC	BANANA	No Resistance	
			KV-0302X146	AHA	SNAP	20 kΩ Resistance	
			KV-0302X147	AHA	CLIP	20 kΩ Resistance	
compatibile: ECG-9320 ECG-9522P (Original part No.: BR-911D)		KV-0302X148	AHA	DIN3.0	20 kΩ Resistance		
		KV-0302X149	AHA	BANANA	20 kΩ Resistance		
		KV-0302X150	IEC	SNAP	20 kΩ Resistance		
		KV-0302X151	IEC	CLIP	20 kΩ Resistance		
		KV-0302X152	IEC	DIN3.0	20 kΩ Resistance		
		KV-0302X153	IEC	BANANA	20 kΩ Resistance		
		PETAS		KV-0302X154	AHA	SNAP	No Resistance
KV-0302X155	AHA			CLIP	No Resistance		
KV-0302X156	AHA			DIN3.0	No Resistance		
KV-0302X157	AHA			BANANA	No Resistance		
KV-0302X158	IEC			SNAP	No Resistance		
KV-0302X159	IEC			CLIP	No Resistance		
KV-0302X160	IEC			DIN3.0	No Resistance		
KV-0302X161	IEC			BANANA	No Resistance		
PHILIPS (HP)		KV-0302X162	AHA	BANANA	20 kΩ Resistance		
		KV-030184	IEC	BANANA	20 kΩ Resistance		
		KV-0302085/1	IEC	CLIP	20 kΩ Resistance		
		KV-0302085	IEC	SNAP	20 kΩ Resistance		
Spacelabs		KV-0302X163	AHA	SNAP	10 kΩ Resistance		
		KV-0302X164	AHA	CLIP	10 kΩ Resistance		
		KV-0302X165	AHA	DIN3.0	10 kΩ Resistance		
		KV-0302X166	AHA	BANANA	10 kΩ Resistance		
		KV-0302070	IEC	SNAP	10 kΩ Resistance		
		KV-0303071	IEC	CLIP	10 kΩ Resistance		
		KV-0302X167	IEC	DIN3.0	10 kΩ Resistance		








LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI

Shanghai Konden 上海光电 ZhonLian JianYi 中联健益	compatibile ShangHai Kohden: ECG-6511, 6151, 6551, 8110P/K, 6353,6403, 6504 JIANYI: XD-106,XD-104B/D		KV-0302X168	IEC	BANANA	10 kΩ Resistance
			KV-0302X169	AHA	SNAP	No Resistance
			KV-0302X170	AHA	CLIP	No Resistance
			KV-0302X171	AHA	DIN3.0	No Resistance
			KV-0302X172	AHA	BANANA	No Resistance
			KV-0302X173	IEC	SNAP	No Resistance
			KV-0302X174	IEC	CLIP	No Resistance
			KV-0302X175	IEC	DIN3.0	No Resistance
			KV-0302X176	IEC	BANANA	No Resistance
			KV-0302X177	IEC	SNAP	No Resistance
Trentina	compatibile: tutti i modelli Excel		KV-0303089/4G	IEC	CLIP	No Resistance
			KV-0302X178	IEC	DIN3.0	No Resistance
			KV-0301089/4F	IEC	BANANA	No Resistance
ZOLL			KV-0301069/2A	IEC	BANANA	10 kΩ Resistance







EKG Multi-link Cable and Lead wires

Compatible Brand	Compatible Models No.	OEM P/N	TH CODE	DESCRIPTION (Name,Length,Plug,STANDARD)	INSTRUMENT CONNECTOR	YOKE
Fukuda			KV-0302X179	12-Lead ECG Trunk Cable, 12pin>>2p*10 Euro Yoke, L=2.2m, AAMI/IEC, 20kΩ Res. -discharge Fit Euro type leadwires		
Hellige GE-Marquette	MicroSmart, MAC 500 / 1100 / 1200 / 1200 ST	223 41808	KV-0302X180	MAC1200 10Ld EKG Trunk Cable,L=2.2m, DB15M>VS 10L Yoke AAMI/IEC, 10kΩ Res		
GE-Medical	compatibile GE Eagle, Solar, Dash , Tram,MAC-Lab Cath Lab System, Datex-Ohmeda S/S FM,Original P/N:(AAMI: E9002ZH, 416035-001) (IEC: E9002ZJ, 416035-002) connettore 11 pin rettangolare		KV-030169/3A	Multi-Link 12-Ld ECG cable L=2.2m, Rectangle 11p>VS 10L Yoke AAMI/IEC, 1kΩ Res.		








LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI

Avionics(Del Mar) Custo-med Del mar Ergoline Esaote Biomedica Fidelity Innomed Kontron Reynolds Schiller Trentina (Elletronica) Von Berg Bionet Welch allyn,	compatibile Avionics(Del Mar): 910/920/930 /940/942/950/960/970/980 , Custo-med: Custo-norm, Ergoline: Autoruler, Autoscript 6/12, Cardiette, EK 3003 / 3012, Ergoline; Esaote Biomedica: EKG P120, Fidelity: LP 3000, Innomed: Heart Mirror(3-Kanal), Kontron: EKG P80, Reynolds: Cardi-track 12, Schiller: AT3, AT5, AT6,AT10, AT60, AT104 PC, AT104, CS6, CS12,CS100, CS200 Trentina: All Excel-Types, Von Berg: Bioset 9000		KV-0302X181	Schiller 10Ld EKG Trunk Cable,L=2.2m DB15M>VS 10L Yoke AAMI/IEC, no Res		
Mindray	compatibile PM-7000/PM-9000Super/PM-9000Express (12 LEAD ECG MODULE) BeneView Oem p/n:0010-30-42721 (AAMI) Oem p/n:0010-30-42722 (IEC)	0010-30-42721	KV-0302X182	Mindray 10Ld EKG Trunk Cable,L=2.2m Round 12p>VS 10L Yoke AAMI/IEC ,1kΩ Res.		
Spacelabs	compatibile TRU-LINK Accessories for 90496 Module	700-0008-01	KV-0302X183	Spacelabs 10Ld EKG Trunk Cable,L=2.2m Round 17p>VS 10L Yoke AAMI/IEC , 1 kΩ Res.		








LEADWIRES

Compatible Brand	Compatible Models No.	OEM P/N	TH CODE	DESCRIPTION (Name,Length,Plug)	STANDARD	PATIENT END TERMINATION
Hellige GE-Maroutte 	compatibile con tutti i MultiLink-Plug System 	384 01816	KV-S03020X012	10 vie Set patient leads, 4 x limb lead (130cm), 6 x chest leads (70cm), no Resistance,VS-2P plug	AHA	SNAP (bottone)
			KV-S03020X013		AHA	BANANA
			KV-S03020X014		IEC	SNAP (bottone)
			KV-S03020X015		IEC	CLIP (pinzetta)
			KV-S03020X016		IEC	BANANA
GE-MEDICAL  	Holter Multi-Link Leadwire Sets – compatibile: GE SEER MC Holter Recorder GE Eagle, Solar, Dash Monitors, Tram, MAC-Lab Cath Lab System, Datex-Ohmeda S/5 FM 		KV-S02020X077	Multi-Link Leadwire Set , 5Ld,L=130cm, AHA	AHA	SNAP (bottone)
			KV-S02020X078	Multi-Link Leadwire Set ,5 vie, Snap, V2-V6, AHA, 51 in. (130cm)	AHA	SNAP (bottone)
			KV-S0203010	Multi-Link Leadwire Set , 5Ld,L=130cm, IEC	IEC	SNAP (bottone)
			KV-S02020X079	Multi-Link Leadwire Set , 5Ld, Snap, C2-C6, IEC, 51 in. (130cm)	IEC	SNAP (bottone)
GE medical 	compatibile GE MAC PC, MAC 6, MAC 12, MAC 15, MAC VU, CASE 12, CASE 15, CASE 16, MAX- Personal, MAX-1, CENTRA Systems with AM4 or AM5 Acquisition Module 		KV-S03020X017	10 vie Set patient leads, 4 x limb lead (130cm), 6 x chest leads (70cm),	AHA	SNAP (bottone)
			KV-S03020X018	no Resistance, L-2P plug	AHA	CLIP (pinzetta)



LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI

				KV-S03020X019		AHA		DIN3.0
				KV-S03020X020		IEC		SNAP (bottone)
				KV-S03020X021		IEC		CLIP (pinzetta)
CONMED JINJIANG	compatibile Cardio Beats Jinjiang Electronic Tech Instrument: LEAD2000, GE electrophysiology machine			KV-S03020X022	10 vie Set patient leads, 4 x limb lead (130cm), 6 x chest leads (70cm), no Resistance, L-2P plug	AHA		SNAP (bottone)
			KV-S03020X023	AHA			DIN3.0	
			KV-S03020X024	IEC			CLIP (pinzetta)	
			KV-S03020X025	IEC			DIN3.0	
			KV-S03020X026	IEC			BANANA	
PHILIPS medical	compatibile Philips Instrument M1700A M1701A M1702A; TRIM II, Trim 1, 2, 3		M1713B	KV-S03020X027		AHA		BANANA
	552		M1716B	KV-S03020X028		IEC		BANANA
Fukuda				KV-S02020X080	Euro style Lead wire Set, 5-Ld, 40in.(102),(RA,LA,RL,LL,V)	AHA		SNAP (bottone)
				KV-S02020X081	Euro style Lead wire Set, 5-Ld, 40in.(102),(V2, V3, V4, V5, V6)	AHA		SNAP (bottone)
				KV-S0203069/2CA	Euro style Lead wire Set, 5-Ld, 40in.(102),(RF,F,C, R,L)	IEC		SNAP (bottone)
				KV-S0202069/2CB	Euro style Lead wire Set, 5-Ld, 40in.(102),(C2, C3, C4, C5, C6)	IEC		SNAP (bottone)

ADAPTER





Compatible Brand	Compatible Models No.	OEM P/N	TH CODE	DESCRIPTION (Name,Length,Plug)	PRODUCTOR GUIDE	PRICE CODE
Hellige GE-Marquette	compatibile GE MAC 1100/1200/1200ST, CardioSmart, CardioSys, MAC 5000, MAC 5000ST, CASE, CASE 8000 Systems with CAM 14 Acquisition Module or with Cable Accepting CAM 14 Leadwires compatibile GE MAC PC, MAC 6, MAC 12, MAC 15, MAC VU, CASE 12,CASE 15, CASE 16, MAXPersonal,MAX-1, CENTRA Systems with AM4 or AMS Acquisition Module	E9006PH		Leadwire,CAM/Universal, 20in.(51cm),6pcs/pk		C4/pk
		E9006PJ		Leadwire,CAM/Universal, 29in.(74cm),6pcs/pk		C4/pk
		E9006PK		Leadwire,CAM/Universal, 40in.(102cm),6pcs/pk		C4/pk
		E9006PL		Leadwire,CAM/Universal, 51in.(130cm),4pcs/pk		C4/pk
		E9006PM		Leadwire,CAM/Universal, 26in.(66cm),6pcs/pk		C4/pk
		E9006PN		Leadwire,CAM/Universal, 36in.(90cm),6pcs/pk		C4/pk
		E9001ZT		Leadwire,AM4/Universal, 24 in.(60cm),6pcs/pk		C4/pk
		E9001ZY		Leadwire,AM4/Universal, 36 in.(90cm),6pcs/pk		C4/pk
		E9001ZZ		Leadwire,AM4/Universal, 42 in.(107cm),6pcs/pk		C4/pk
		E9002A	KV-S030201AM4-5	Leadwire,AM4/Universal, 52 in.(132cm),4pcs/pk		143

LISTINO CAVI ECG COMPATIBILI



PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
il presente listino può essere soggetto a variazioni senza preavviso

EKG Multi-Link Cable and Leadwires











EKG Multi-Link Cable

Compatible Brand	Compatible Models No.	OEM P/N	TH CODE	DESCRIPTION (Name,Length,Plug,STANDARD)	INSTRUMENT CONNECTOR	YOKE
MEDTRONIC PHYSIO CONTRO	compatibile Medtronic Physio Control: Lifepak 12		KV-0302X184	12-Lead ECG Trunk Cable, 12pin->LL4+7pin Yoke, L=2.0m, AAMI/IEC, 10kΩ Res. Fit LL 4LD leadwires and VR02456A, VR02456I		
Spacelabs	compatibile 90496 Ultraview Module		KV-0302X185	12-Lead ECG Trunk Cable, 17pin->LL4+7pin Yoke, L=2.0m, AAMI/IEC, 1kΩ Res Fit LL 4LD leadwires and VR02456A, VR02456I		

EKG Multi-Link Lead set

Compatible Brand	Compatible Models No.	OEM P/N	TH CODE	DESCRIPTION (Name,Length,Plug,STANDARD)	DIAGRAM
EKG Multi-Link Cable	compatibile Medtronic Physio Control:Lifepak 12, Zoll medical:R Series Monitor Def.		KV-S02020X82	AAMI 6 vie ECG leadwires, 7pin->6*Snap, L=1.15m,	
			KV-S02020X83	IEC 6 vie ECG leadwires, 7pin->6*Snap, L=1.15m,	
EKG Multi-Link Cable	compatibile Medtronic Physio Control:Lifepak 12, Zoll medical:R Series Monitor Def.		KV-S02010X61	AAMI 4 vie ECG leadwires, LL 2pin >> 4*Snap, L=1.0m,	
			KV-S02010X62	IEC 4 vie ECG leadwires, LL 2pin >> 4*Snap, L=1.0m,	

















Holter Trunk Cables

Compatible Brand	Compatible Models No.	OEM P/N	TH CODE	DESCRIPTION (Name,Length,Plug,STANDARD)	INSTRUMENT CONNECTOR	YOKE
Biomedical Instruments (BI)	compatibile Biomedical Instruments (BI): BI 9000 Brentwood : BH-3000		KV-H02010XX1	Din 7-Ld Patient Cable, 3CH, 45 cm , AHA, FGG.0B 7pin->Din 7j, no resistance fit for Din style LeadWires		
	compatibile BI 9800		KV-H02010XX2	Din 10 vie Patient Cable, 45 cm , AHA, Reader-14pin->Din 10j, no resistance fit for Din style LeadWires		
Biox	compatibile Biox: CB-2300-A Biomedical Instruments (BI): BI9000TL		KV-H02010XX3	Din 7-Ld Patient Cable, 3CH, 45 cm , AHA, Reader-14pin->Din 7j, no resistance fit for Din style LeadWires		
Brentwood	compatibile Biomedical Instruments (BI): BI 9000 Brentwood : BH-3000		KV-H02010XX4	Din 7-Ld Patient Cable, 3CH, 45 cm , AHA, FGG.0B 7pin->Din 7j, no resistance fit for Din style LeadWires		
Burdick	compatibile 92512 Holter	43238	KV-H02010XX5	Din 5 vie Patient Cable, 2CH, 80 cm , AHA, PCB 7pin->Din 5j, no resistance fit for Din style LeadWires		



LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI

	compatibile 92513 Holter	175-1546-00	KV-H02010XX6	Din 5 vie Patient Cable, 3CH, 80 cm , AHA, PCB 7pin>Din 5J, no resistance fit for Din style LeadWires		
Landcom			KV-H02010XX7	Din 7-Ld Patient Cable, 3 CH, 45 cm , AHA, FGG.1B 7pin>Din 7J, no resistance fit for Din style LeadWires		
			KV-H02010XX8	Din 10 vie Patient Cable, 45 cm , AHA, FGG.1B 10pin>Din 10J, no resistance fit for Din style LeadWires		
Medilog	compatibile Medilog: FD5		KV-H02010XX9	Din 7-Ld Patient Cable, 3CH, 45 cm , AHA, Readel-14pin>Din 7J, no resistance fit for Din style LeadWires		
Ge-Marquette Rozinn GeTeMED	compatibile GE marquette: 8000, 8200, 8500, SEER Holter Recorder model: 151	E9002TE	KV-0101MQ01	Snap-Lock Holter Cable, 3 CH, 45 cm , AHA, Snap-lock-7pin>Din 7J, no resistance fit for Din style LeadWires		
GE-Marquette	compatibile GE SEER MC Holter Recorder	E9002TK	KV-H02010X10	Din 5 vie Patient Cable, 2 CH, 80 cm , AHA, rectangle 12pin>Din 5J, no resistance fit for Din style LeadWires		
		E9002TL	KV-H02010X11	Din 7-Ld Patient Cable, 3 CH, 80 cm , AHA, rectangle 12pin>Din 7J, no resistance fit for Din style LeadWires		
Rozinn		E9002TP	KV-H02010X12	Din 10 vie Patient Cable, 10 vie, 80 cm , AHA, rectangle 12pin>Din 10J, no resistance fit for Din style LeadWires		












One-piece ECG patient cable with lead wires

Compatible Brand	Compatible Models No.	OEM P/N	TH CODE	DESCRIPTION (Name,Length,Plug,STANDARD)	STANDARD	CONNECTOR
Biomedical Instruments (BI)	B1900TL 		KV-H02010X13	One-piece holter ECG cable with leadwires, 7-Ld Plug: Radel 14pin	AAMI	SNAP (bottone)
Choice Telemetry Monitoring System MMED600RD-300E			KV-H02010X14	3Ld lead wires, 40 in. (102 cm) , 3 vie PLUG: READEL 6P	IEC	SNAP (bottone)
			KV-H02010X15	3Ld lead wires, 40 in. (102 cm) , 3 vie PLUG: READEL 5P	IEC	CLIP (pinzetta)

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso





LISTINO PREZZI RICAMBI COMPATIBILI

Creative PC-6000 Telemetry Monitoring			KV-H02010X16	3Ld lead wires, 40 in. (102 cm), 3 vie PLUG: Din2.0*3pin	AHA		SNAP (bottone)
			KV-H02010X17	3Ld lead wires, 40 in. (102 cm), 3 vie PLUG: stereo 2.5 plug	AHA		SNAP (bottone)
Creative PC-80 Easy ECG Monitor			KV-H02010X18	3Ld lead wires, 40 in. (102 cm), 3 vie PLUG: stereo 2.5 plug	AHA		SNAP (bottone)
			KV-H02010X19	3Ld lead wires, 40 in. (102 cm), 3 vie PLUG: stereo 2.5 plug	AHA		SNAP (bottone)
Fukuda	LX-5630 		KV-H02010X20	One-piece holter ECG cable with leadwires, 4 vie Plug: 6 pin	AAMI		CLIP (pinzetta)
Meigaoyi MGY-M1 Landcom Visual Stethoscope			KV-H02010X21	3Ld lead wires, 40 in. (102 cm), 3 vie PLUG: stereo 3.5 plug	IEC		SNAP (bottone)
			KV-H02010X22	3Ld lead wires, 40 in. (102 cm), 3 vie PLUG: stereo 3.5 plug	IEC		CLIP (pinzetta)
Meigaoyi ECGLAB TDA3.0			KV-H02010X23	5 vie lead wires, 40 in. (102 cm), 5 vie PLUG: stereo 2.5+3.5 plug	AHA		SNAP (bottone)
Mortara			KV-0401M1	10 vie lead wires, 40 in. (102 cm)	IEC		SNAP (bottone)
GE-Marquette Rozinn			KV-H02010X24	3 vie lead wires, 40 in.	AHA		SNAP (bottone)
GE-SEER Light		2008594-001	KV-H02010X25	SEER Light Holter Patient Cable/Leadwires, 5 vie	AHA		SNAP (bottone)
		2008594-002	KV-H02010X26	SEER Light Holter Patient Cable/Leadwires, 7-Ld	AHA		SNAP (bottone)
		2008594-004	KV-H02010X27	SEER Light Holter Patient Cable/Leadwires, 7-Ld	IEC		SNAP (bottone)
Mindray TMS-6016 Telemetry monitoring system			KV-H02010X28	5Ld lead wires, 40 in. (102 cm), 5 vie PLUG: 1.5mm*5pin	AHA		SNAP (bottone)
Mbelec CD2000-1800			KV-H02010X29	3Ld lead wires, 40 in. (102 cm), 3 vie PLUG: 1.5mm*5pin	AHA		SNAP (bottone)
NEC		YCE-205	KV-H02010X30	YCE205 lead wires, 30 in. (76 cm), 3 vie	IEC		SNAP (bottone)
			KV-H02010X31	YCE205 lead wires, 30 in. (76 cm), 3 vie	IEC		CLIP (pinzetta)



LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI

		47502	KV-H02010X32	47502 lead wires, 30 in. (76 cm), , 3 vie	IEC		SNAP (bottone)
Nihon kohden	 124	BR-546S	KV-H02010X33	3Ld lead wires, 40 in. (102 cm) , 3 vie PLUG: Din2.0*3pin	IEC		SNAP (bottone)
			KV-H02010X34	3Ld lead wires, 40 in. (102 cm) , 3 vie PLUG: Din2.0*3pin	IEC		CLIP (pinzetta)
PI			KV-H02010X35	4 vie lead wires,40 in. (102 cm) , PLUG: HRS-4P	IEC		SNAP (bottone)
			KV-H02010X36	6 vie lead wires,40 in. (102 cm) , PLUG: HRS-6P	IEC		SNAP (bottone)
Philips			KV-H02010X37	Telemetry monitoring leads ,40 in. (102 cm) ,	IEC		SNAP (bottone)
			KV-H02010X38	Telemetry monitoring leads ,40 in. (102 cm) ,	IEC		SNAP (bottone)
Din 1.5 Style Normal	compatibile con tutti i Din-Plug System , Yoke type: Din style,Code: C		KV-H02010X39	Din lead wires , 18 in. (45 cm) , 3 vie	AHA		SNAP (bottone)
			KV-H02010X40	Din lead wires , 18 in. (45 cm) , 3 vie	AHA		CLIP (pinzetta)
			KV-H02010X41	Din lead wires , 18 in. (45 cm) , 3 vie	IEC		SNAP (bottone)
			KV-H02010X42	Din lead wires , 18 in. (45 cm) , 3 vie	IEC		CLIP (pinzetta)
			KV-H02010X43	Din lead wires , 18 in. (45 cm) , 5 vie	AHA		SNAP (bottone)
			KV-H02010X44	Din lead wires , 18 in. (45 cm) , 5 vie	AHA		CLIP (pinzetta)
			KV-H02010X45	Din lead wires , 18 in. (45 cm) , 5 vie	IEC		SNAP (bottone)
			KV-H02010X46	Din lead wires , 18 in. (45 cm) , 5 vie	IEC		CLIP (pinzetta)






LEADWIRES (TERMINALI)

Compatible Brand	Compatible Models No.	OEM P/N	TH CODE	DESCRIPTION (Name,Length,Plug)	STANDARD	PATIENT END TERMINATION
Din 1.5 Style Normal	compatibile con tutti i DIN-plug system  141  EC024S3A		KV-H02010X47	Din lead wires,24 in. (61 cm) , 3 vie	AHA	SNAP (bottone)
			KV-H02010X48	Din lead wires,24 in. (61 cm) , 3 vie	AHA	CLIP (pinzetta)
			KV-H02010X49	Din lead wires,24 in. (61 cm) , 3 vie	IEC	SNAP (bottone)
			KV-H02010X50	Din lead wires,24 in. (61 cm) , 3 vie	IEC	CLIP (pinzetta)
			KV-H02010X51	Din lead wires,24 in. (61 cm) , 5 vie	AHA	SNAP (bottone)
			KV-H02010X52	Din lead wires,24 in. (61 cm) , 5 vie	AHA	CLIP (pinzetta)
			KV-H02010X53	Din lead wires,24 in. (61 cm) , 5 vie	IEC	SNAP (bottone)
			KV-H02010X54	Din lead wires,24 in. (61 cm) , 5 vie	IEC	CLIP (pinzetta)
			KV-H02010X55	Din lead wires,40 in. (102 cm) , 3 vie	AHA	SNAP (bottone)
			KV-H02010X56	Din lead wires,40 in. (102 cm) , 3 vie	AHA	CLIP (pinzetta)
			KV-H02010X57	Din lead wires,40 in. (102 cm) , 3 vie	IEC	SNAP (bottone)
			KV-H02010X58	Din lead wires,40 in. (102 cm) , 3 vie	IEC	CLIP (pinzetta)
			KV-H02010X59	Din lead wires,40 in. (102 cm) , 5 vie	AHA	SNAP (bottone)
			KV-H02010X60	Din lead wires,40 in. (102 cm) , 5 vie	AHA	CLIP (pinzetta)
			KV-H02010X61	Din lead wires,40 in. (102 cm) , 5 vie	IEC	SNAP (bottone)
			KV-H02010X62	Din lead wires,40 in. (102 cm) , 5 vie	IEC	CLIP (pinzetta)

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI

		KV-H02010X63	Din lead wires,18 in. (45 cm) , 7-Ld	AHA	SNAP (bottone)
		KV-H02010X64	Din lead wires,18 in. (45 cm) , 7-Ld	AHA	CLIP (pinzetta)
		KV-H02010X65	Din lead wires,18 in. (45 cm) , 7-Ld	IEC	SNAP (bottone)
		KV-H02010X66	Din lead wires,18 in. (45 cm) , 7-Ld	IEC	CLIP (pinzetta)
		KV-H02010X67	Din lead wires,24 in. (61 cm) , 7-Ld	AHA	SNAP (bottone)
		KV-H02010X68	Din lead wires,24 in. (61 cm) , 7-Ld	AHA	CLIP (pinzetta)
		KV-H02010X69	Din lead wires,24 in. (61 cm) , 7-Ld	IEC	SNAP (bottone)
		KV-H02010X70	Din lead wires,24 in. (61 cm) , 7-Ld	IEC	CLIP (pinzetta)
		KV-H02010X71	Din lead wires,40 in. (102 cm) , 7-Ld	AHA	SNAP (bottone)
		KV-H02010X72	Din lead wires,40 in. (102 cm) , 7-Ld	AHA	CLIP (pinzetta)
		KV-H02010X73	Din lead wires,40 in. (102 cm) , 7-Ld	IEC	SNAP (bottone)
		KV-H02010X74	Din lead wires,40 in. (102 cm) , 7-Ld	IEC	CLIP (pinzetta)
		KV-H02010X75	Din lead wires,18 in. (45 cm) , 10 vie	AHA	SNAP (bottone)
		KV-H02010X76	Din lead wires,18 in. (45 cm) , 10 vie	AHA	CLIP (pinzetta)
		KV-H02010X77	Din lead wires,18 in. (45 cm) , 10 vie	IEC	SNAP (bottone)
		KV-H02010X78	Din lead wires,18 in. (45 cm) , 10 vie	IEC	CLIP (pinzetta)
		KV-H02010X79	Din lead wires,24 in. (61 cm) , 10 vie	AHA	SNAP (bottone)
		KV-H02010X80	Din lead wires,24 in. (61 cm) , 10 vie	AHA	CLIP (pinzetta)
		KV-H02010X81	Din lead wires,24 in. (61 cm) , 10 vie	IEC	CLIP (pinzetta)
		KV-H02010X82	Din lead wires,40 in.(102 cm) , 10 vie	AHA	SNAP (bottone)
	KV-H02010X83	Din lead wires,40 in.(102 cm) , 10 vie	AHA	CLIP (pinzetta)	
	KV-H02010X84	Din lead wires,40 in.(102 cm) , 10 vie	IEC	CLIP (pinzetta)	
Contec, Meigaoyi	 <p>133 fit 1.35~1.6mm pin</p>	KV-H02010X85	Din lead wires,40 in. (102 cm) , 7-Ld	AHA	SNAP (bottone)
		KV-H02010X86	Din lead wires,40 in.(102 cm) , 10 vie	AHA	SNAP (bottone)
Din 2.0 Jack Style Biomedical system: DH 073	compatibile con tutti i DIN 2.0 jack system yoke type: DIN 2.0, Code: S  <p>32 φ2.0 mm Jack</p>	KV-H02010X87	Din 2.0 lead wires,24 in. (61 cm) , 7-Ld	AHA	SNAP (bottone)
		KV-H02010X88	Din 2.0 lead wires,24 in. (61 cm) , 7-Ld	AHA	CLIP (pinzetta)
		KV-H02010X89	Din 2.0 lead wires,24 in. (61 cm) , 7-Ld	AHA	CLIP (pinzetta)
Din 2.0 Plug Style Ela Medical	compatibile con tutti i DIN 2.0 plug system yoke type: DIN 2.0 plug style, Code: T  <p>134</p>	KV-H02010X90	Din 2.0 plug lead wires,40 in. (102 cm) , 3 vie	AHA	SNAP (bottone)
		KV-H02010X91	Din 2.0 plug lead wires,40 in. (102 cm) , 3 vie	IEC	CLIP (pinzetta)
		KV-H02010X92	Din 2.0 plug lead wires,40 in. (102 cm) , 5 vie	IEC	CLIP (pinzetta)
LL1 Plug Style Medilog: FD4	compatibile con tutti i LL1-plug system yoke type: DIN style, Code: U  <p>337</p>	KV-H02010X93	18 in. (45 cm) , 3 vie	AHA	SNAP (bottone)
		KV-H02010X94	18 in. (45 cm) , 3 vie	AHA	CLIP (pinzetta)
		KV-H02010X95	18 in. (45 cm) , 3 vie	IEC	SNAP (bottone)
		KV-H02010X96	18 in. (45 cm) , 3 vie	IEC	CLIP (pinzetta)
		KV-H02010X97	18 in. (45 cm) , 5 vie	AHA	SNAP (bottone)
		KV-H02010X98	18 in. (45 cm) , 5 vie	AHA	CLIP (pinzetta)
		KV-H02010X99	18 in. (45 cm) , 5 vie	IEC	SNAP (bottone)
		KV-H02010X100	18 in. (45 cm) , 5 vie	IEC	CLIP (pinzetta)
		KV-H02010X101	24 in. (61 cm) , 3 vie	AHA	SNAP (bottone)
		KV-H02010X102	24 in. (61 cm) , 3 vie	AHA	CLIP (pinzetta)
		KV-H02010X103	24 in. (61 cm) , 3 vie	IEC	SNAP (bottone)
		KV-H02010X104	24 in. (61 cm) , 3 vie	IEC	CLIP (pinzetta)
		KV-H02010X105	24 in. (61 cm) , 5 vie	AHA	SNAP (bottone)
		KV-H02010X106	24 in. (61 cm) , 5 vie	AHA	CLIP (pinzetta)
		KV-H02010X107	24 in. (61 cm) , 5 vie	IEC	SNAP (bottone)
		KV-H02010X108	24 in. (61 cm) , 5 vie	IEC	CLIP (pinzetta)
		KV-H02010X109	40 in. (102 cm) , 3 vie	AHA	SNAP (bottone)
		KV-H02010X110	40 in. (102 cm) , 3 vie	AHA	CLIP (pinzetta)
		KV-H02010X111	40 in. (102 cm) , 3 vie	IEC	SNAP (bottone)
		KV-H02010X112	40 in. (102 cm) , 3 vie	IEC	CLIP (pinzetta)
		KV-H02010X113	40 in. (102 cm) , 5 vie	AHA	SNAP (bottone)

LISTINO CAVI ECG COMPATIBILI

PS I codici con "X" sono da ritenersi provvisori, richiedere codice presso nostri uffici
 il presente listino può essere soggetto a variazioni senza preavviso

LISTINO PREZZI RICAMBI COMPATIBILI



KV-H02010X114	40 in. (102 cm) , 5 vie	AHA	CLIP (pinzetta)
KV-H02010X115	40 in. (102 cm) , 5 vie	IEC	SNAP (bottone)
KV-H02010X116	40 in. (102 cm) , 5 vie	IEC	CLIP (pinzetta)
KV-H02010X117	18 in. (45 cm) , 7-Ld	AHA	SNAP (bottone)
KV-H02010X118	18 in. (45 cm) , 7-Ld	AHA	CLIP (pinzetta)
KV-H02010X119	18 in. (45 cm) , 7-Ld	IEC	SNAP (bottone)
KV-H02010X120	18 in. (45 cm) , 7-Ld	IEC	CLIP (pinzetta)
KV-H02010X121		AHA	CLIP (pinzetta)
KV-H02010X122		IEC	SNAP (bottone)
KV-H02010X123		IEC	CLIP (pinzetta)
KV-H02010X124	40 in. (102 cm) , 7-Ld	AHA	SNAP (bottone)
KV-H02010X125	40 in. (102 cm) , 7-Ld	AHA	CLIP (pinzetta)
KV-H02010X126	40 in. (102 cm) , 7-Ld	IEC	SNAP (bottone)
KV-H02010X127	40 in. (102 cm) , 7-Ld	IEC	CLIP (pinzetta)

Condizioni di fornitura:

- Ai prezzi sopra riportati va applicata un'aliquota IVA pari al 21%.
- Tempo di consegna: entro 30 giorni solari dalla data di ricezione del Vostro ordine.
- Spedizione: a Nostro carico
- Modalità di pagamento: Rim. Dir. 60gg D.F. oltre a tale scadenza verranno applicati interessi come legislazione vigente
- Validità offerta: fino al 31/12/2012.
- Minimo fatturabile: €. 150,00 + IVA 21%
- In caso di ordine, attenersi ai confezionamenti eventualmente indicati, **citare sempre il nostro riferimento ed i nostri codici** e specificare i vostri dati (sede fiscale/legale, coordinate bancarie, P.I., codice fiscale) nonché il luogo di consegna, gli orari di ricevimento della merce ed altre informazioni utili.
- La presente offerta non vincola in alcun modo Tecnohealth s.r.l. la quale si riserva il diritto di non accettare ordini ad essa riferiti.
- **NOTA BENE:** i codici con la **X** sono da ritenere provvisori (es. KV-M02010X001)